# M Bike motorise

Réf : 86075

F

GB

D

E

I

NL

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### Sécurité

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## Important: Read these instructions before using the equipment!

### Safety

- Read this user manual and keep it in a safe place. Use this product exclusively in the manner indicated.
- This equipment must be assembled and used by adults.
- The equipment complies with the EN-957 standard, classification H for use in the home. It must not be used in a sports center or any other public, associative or rented space.
- For totally safe use, a stable, level surface is required. Protect your floor covering with a mat. Do not use the equipment in damp areas (swimming pool, sauna, etc.).
- Do not allow children to play with this equipment. CARE declines all responsibility for any injuries they may incur. Do not allow children close to this equipment when you are training.
- CARE declines all responsibility if technical modifications are made to any of its products by the user.
- Maximum weight of user: 80 kg.
- **Before you start training, it is essential that you consult your doctor to determine the intensity level of the program you use.**
- Excessive or badly programmed training can damage your health.
- Keep your back straight when exercising.
- It is strongly recommended that you wear suitable clothing and footwear.

### Maintenance

- Regularly check that elements fastened with nuts and bolts are correctly tightened.
- To maintain the level of safety, the equipment must be regularly checked. It is vital that any faulty parts are replaced and that it is not used until completely repaired. Remember regularly to grease moving parts.
- As sweat is very corrosive, do not allow it to come into contact with the enameled or chromed parts of the equipment, particularly the computer. Immediately wipe the equipment after training. The enameled parts can be cleaned using a damp sponge. All aggressive or corrosive products should be avoided.

**Guarantee:** the chassis is guaranteed for 5 years. Moving parts are guaranteed for 2 years. The guarantee applies to normal use by a private individual in his home. Validate your product warranty on line [www.carefitness.com](http://www.carefitness.com)

**Recycling :** The « crossed out dustbin » sign means that this product and its batteries cannot be thrown out with domestic waste. They should be treated apart. When you have finished with them, drop them at an authorised collection point so they can be recycled. This gesture will go towards protecting the environment and your health.

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E mail : [sav@carefitness.com](mailto:sav@carefitness.com)

## ASSEMBLY INSTRUCTIONS.-

Take the unit out of the box and place it on the floor as shown in Fig.1.

### The assistance of a second person is advisable for the assembly work.

Make sure that all of the pieces are there:

(A) Main body;(C-L) left pedal(C-R) Right pedal;(B)Stabiliser bar;(B1) Rear stabiliser bar;(D) Screw M-4X25; (E) Grower washer M-4; Pedal combination spanner.

### 1.-FITTING THE STABILISERS.-

Take the stabiliser (B) marked with the letter(A) and place it at the front of the unit marked with the letter (A), as shown in Fig.2, insert the M-4X25 screws(D).

Next, take the rear stabiliser (B1) marked with the letter (B), Fig.3, and go through the same fitting procedure as with the front stabiliser.

### 2. ATTACHING THE PEDALS.-


Right and left refer to the position that the user adopts when to do the exercises. Insert the right-hand pedal (C-R), marked with the letter (R), onto the right-hand crank, also marked with an (R). Next, Fig.4.

Insert the left-hand pedal (C-L), marked with the letter (L), onto the left-hand crank, also marked with an (L), Fig.4.

### 3.MAINS CONNECTION.-

Connect the unit to a plug provided of a ground (220v). Put the switch (L) IN POSITION [1] Fig.6.

## REMOTE CONTROL

**Display indicator:** Speed, Time, Distance, (Distance traveled km ) Calorie, Scan (The random sample data display speed, time, distance, calories.)

- 1: **Power**
- 2.: Time-Time selection button.(1-15min).
- 3: Speed - : Select key for increasing speed (12 speed levels).
- 4: Speed: Select key for increasing speed (12 speed levels)
- 5: Forward/reverse- change pedal direction key
- 6: Mode-select key for parameter display speed, time , distance, calories.


## EXERCICES

### Position I

#### Seated . -

##### lower extremities


1. Use a stable and sturdy chair.
2. Be sure to place the Mini-Bike at a comfortable distance from the chair (so as not to overstretch your legs).
3. Make sure that your lower back (lumbar region) is supported when using the Mini-Bike. We recommend placing a small cushion at the base of the lumbar region.
4. Sit and place your feet on the pedals.
5. Take the control unit off of its support and press the button (1), the pedals will begin to rotate.
6. To stop the unit press (1).

### Position II

#### Lying . -

1. When lying down make sure that your lower back (lumbar region) is supported when using the Mini-Bike. We recommend placing a small cushion below the lumbar region and under your head.
2. Lie with both feet on the pedals (you can use the support straps if you wish).
3. Place your arms at your side and follow the steps in Position I .


4. When you finish the session, you can continue lying in the same position for a few moments and then get up slowly.

**Important Note if you are using the appliance on the bed make sure that the Mini-Bike is stable. Place it on top of a board or stable support so that it does not rest directly on the mattress or bedclothes.**

### Position III

#### Seated. - upper body

1. Be sure to place the Mini-Bike on a stable surface and at a comfortable distance from the chair (so as not to overstretch your arms).
2. Place the Mini-Bike on a sturdy table (for example) and sit on a stable chair. Position the Mini-Bike at chest height with the support for the control unit facing you.
3. Follow the steps for Position I to begin the session.
4. When you finish the session, stay in the same position for a few moments and then get up slowly.

### After Sales Service

Email: [service-uk@carefitness.com](mailto:service-uk@carefitness.com)

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