SANITAS

SFT 65

D	MULTIFUNKTIONS-THERMOMETER Gebrauchsanweisung	2-17
GB	MULTI-FUNCTION THERMOMETER Instructions for use	18-31
F	THERMOMÈTRE MULTIFONCTION Mode d'emploi	32-45
I	TERMOMETRO MULTIFUNZIONE Istruzioni per l'uso	46-59
RUS	МНОГОФУНКЦИОНАЛЬНЫЙ ТЕРМОМЕТР Инструкция по применению	60-75
	ELECTROMAGNETIC COMPATIBILITY INFORMATION	76-79

Service-Adresse:
Hans Dinslage GmbH,
Riedlinger Str. 28,
88524 Uttenweiler, GERMANY
Tel.: +49 (0) 73 74 / 91 57 66,
E-mail: service@sanitas-online.de

Deutsch

Inhalt

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English

Contents

1	Signs and symbols	18	7	Measuring	24
2	Tips	19	8	Replacing the battery.....	27
3	Device description	21	9	Storage and maintenance.....	28
4	Getting started	22	10	Disposal	29
5	Setting the device	22	11	Technical data	29
6	What you should know about using this thermometer	23	12	Trouble-shooting	31

Included in delivery

- Thermometer
- 2 x 1.5V AAA batteries, LR03
- These Instructions for Use

Please read through these Instructions for Use carefully and retain them for future reference; make them available to other users and observe the information they contain.

1 Signs and symbols

The following symbols appear in these instructions.

	Warning	Warning instruction indicating a risk of injury or damage to health.
	Important	Safety note indicating possible damage to the unit/accessory.
	Note	Note on important information.

2 Tips

Safety tips

Warning:

- To measure temperature, insert the sensor tip of the thermometer carefully into the ear.
- Use of the thermometer on different persons can be inappropriate in the event of certain acute infectious diseases because of the possible spread of germs despite cleaning and disinfection. If you have any doubts, please consult your doctor.

Notes on handling batteries

Warning:

- Batteries can be extremely dangerous if swallowed. Keep batteries out of the reach of small children. Should a battery be swallowed, seek medical assistance immediately.
- Batteries must not be recharged or reactivated by any other means, taken apart, thrown into an open fire or short circuited.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable legal regulations. Do not dispose of batteries with the normal household waste.

Important

- Replace weak batteries before they discharge completely.
- Leaking batteries can damage the device. If you do not use the device for an extended period of time, remove the batteries from the battery compartment.
- Leaking batteries may damage the device. If you do not intend to use the device for longer periods, remove the batteries from the battery compartment.



General advice

- Please read through these Instructions for Use carefully and retain them for future reference; make them available to other users and observe the information they contain.
 - This thermometer is a sensitive electronic device. Handle it with care and do not subject it to any mechanical impacts.
 - Do not expose the device to direct sunlight.
 - Check before each use that the lens is intact. If it is damaged, please contact your retailer or the service address.
 - The thermometer is NOT waterproof. For this reason, direct contact with water or other liquids should be avoided.
 - This device is not intended for commercial or clinical use.
 - The SFT 65 thermometer is only designed for the measuring area on the human body stated in the instructions for use.
 - The device is only intended for the purpose stated in these Instructions for Use.
 - After each use, clean the measuring tip with a soft cloth moistened with disinfectant.
 - The thermometer has been designed for practical use, but is not a substitute for a visit to the doctor.
 - Before each repair query, first check the batteries and replace them if necessary.
 - Repairs must only be carried out by authorised service centres, otherwise the guarantee is invalidated.
 - This device complies with the EU Directive 93/42/EC concerning medical products, the Medizinproduktegesetz (German Medical Devices Act), the ASTM (American Society for Testing and Materials) E 1965-98 and the European Standard EN 12470-5: Clinical thermometers – Part 5: Requirements for infrared ear thermometers (with maximum device) and the European Standard EN60601-1-2, and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Service address or found at the end of the instructions for use.
 - If you should have any questions about using our devices, please contact either your retailer or Customer Services.
-

3 Device description


Overview

- 1 Forehead/cover cap
- 2 Forehead thermometer mode button
- 3 Ear thermometer mode/On button /
- 4 LED (red) (high temperature alarm)
- 5 LED (green)

- 11 Date
- 12 Time
- 13 Measurement mode ("Ear", "Forehead", "Object")
- 14 Battery warning symbol
- 15 Measurement display
- 16 "Memory" symbol and storage space number
- 17 "Measurement in progress" symbol
- 18 Unit of measurement ($^{\circ}\text{C}$ or $^{\circ}\text{F}$)

- 6 LCD
- 7 Sensor tip with lens (cover cap removed)
- 8 Memory button
- 9 Set button
- 10 Battery compartment lid


4 Getting started

The batteries are already inserted in all new devices. Before first use, remove the protruding plastic insulating strip from the battery compartment. The thermometer will then switch on automatically.

5 Setting the device

Set the unit in which the temperature will be displayed, the date and the time

Press the **I** button for 1 second to switch the thermometer on. The device performs a brief self test: all elements of the display are shown for approx. 1 second. Following a successful self test, the device emits two short beeps. In Standby mode, the date and time are shown and “- -” is displayed.

When the thermometer is used for the first time, and each time the batteries are changed, the date and time are displayed in the format **I- I 00:00** during the basic set-up (after the self test). Now set the unit in which the temperature will be displayed, the date and the time. To set the device, proceed in stages: unit – year – month – day – hours – minutes.

- Press the **S** button for approx. 3 second to access the Set mode.
The setting to be changed will flash in turn.
- Press the **M** button again to save the setting or
- press the **S** button to change the setting.
- Then press the **M** button to save the setting.

After this process, the thermometer shows the word **OFF** in the display and switches off automatically.

6 What you should know about using this thermometer

There are different thermometers for taking temperature measurements in different parts of the body:

- Ear/forehead thermometer (this thermometer, for measuring in the ear or on the forehead area)
- Rod thermometer (for rectal (in the back passage), axillary (under the arm) or oral measurement (in the mouth))

Taking the temperature gives a measurement that provides information about a person's current body temperature. If you are uncertain about interpreting the results or if the value is abnormal, please consult your doctor. This also applies in the case of slight temperature changes if there are other symptoms of illness such as agitation, severe sweating, flushed skin, fast pulse rate, tendency to collapse etc.

Temperatures measured with different thermometers cannot be compared with one another. You should therefore tell your doctor (or bear in mind if you are diagnosing yourself) what thermometer you used to take your temperature and in what part of the body.

The temperature of a healthy person is affected by various factors: the person's individual metabolism, their age (body temperature is higher in babies and toddlers and decreases with age). Greater temperature fluctuations occur faster and more often in children, e.g. due to growth spurts), their clothing, the ambient temperature, the time of day (body temperature is lower in the morning and increases throughout the day towards evening), the preceding physical and, to a lesser extent, mental activity.

The temperature varies depending on the part of the body where the measurement is taken. The difference can be between 0.4 °F (0.2 °C) and 1.8 °F (1 °C) for a healthy body. The normal temperature ranges are as follows:

- On the forehead: 96.4 °F – 99.7 °F (35.8 °C – 37.6 °C), measured with a forehead thermometer
- In the ear: 96.8 °F – 100 °F (36.0 °C – 37.8 °C), measured with an ear thermometer
- In the back passage: 97.3 °F – 100 °F (36.3 °C – 37.8 °C), measured with a conventional thermometer
- In the mouth: 96.8 °F – 99.3 °F (36.0 °C – 37.4 °C), measured with a conventional thermometer

In order to track temperature change, always take the measurements in the same part of the body.

7 Measuring



Before each use, please check that the lens is intact. If it is damaged, please contact your retailer or the service address. Bear in mind that the thermometer needs to have been in the room in which the measurement is taken for at least 30 minutes before use.

Measuring body temperature in the ear

- Some people produce different readings in their left and right ear. In order to record temperature changes, always measure a person's temperature in the same ear.
- The ear thermometer may be used by children only under adult supervision. Measurement is usually possible over the age of 6 months. In infants under 6 months, the ear canal is still very narrow so the temperature of the eardrum often cannot be recorded and the result displayed is often too low.
- The measurement must not be taken in an ear affected by inflammatory diseases (e.g. discharging pus or secretion), after possible ear injuries (e.g. eardrum damage) or in the healing period after operative procedures. In all of these cases, please consult your doctor.
- Use of the thermometer on different persons can be inappropriate in the event of certain acute infectious diseases because of the possible spread of germs despite cleaning and disinfection. If you have any doubts, please consult your doctor.
- This thermometer may only be used without a disposable protective cover.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.
- As ear wax can affect the measurement, you should clean the ear before measuring if necessary.

- Press the **I** button for 1 second to switch the thermometer on. Following a successful self test, the device emits two short beeps.
- Remove the cap by gently pushing it up (1) and then pulling it off forwards (2).


Make sure that the sensor tip and also the ear canal are clean.

As the ear canal is slightly curved, you have to pull the ear slightly up and backwards before inserting the sensor tip. This is important so that the sensor tip can be pointed directly at the eardrum.

- Insert the sensor tip carefully and press the **I** button for 1 second.
- Release the **I** button. The end of the measuring time is signalled with a short beep and the measured value appears on the display.

If the measured value is within the normal range (< 100.4 °F/38 °C), the green LED lights up for 3 seconds. If the measured value is higher (\geq 100.4 °F/38 °C, high temperature), the red LED is illuminated.

Measuring body temperature on the forehead

Please note that the forehead/temples must be free from sweat and cosmetics and that taking vasoconstrictive medication and skin irritations can distort the result when measuring temperature on the forehead.

- Press the **I** button for 1 second to switch the thermometer on. Following a successful self test, the device emits two short beeps.
- Place the measuring head with the forehead/cover cap fitted on the temples, hold down the **Stirn** button and move the thermometer smoothly over the forehead to the other temple.
- Then release the **Stirn** button. The end of the measuring time is signalled with a short beep and the measured value appears on the display.

If the measured value is within the normal range (< 100.4 °F/38 °C), the green LED lights up for 3 seconds. If the measured value is higher (≥ 100.4 °F/38 °C, high temperature), the red LED is illuminated.

Measuring surface temperature

- Press the **I** button for 1 second to switch the thermometer on. Following a successful self test, the device emits two short beeps.
- Then press the **M** and **I** buttons simultaneously for 3 seconds to switch to the Object temperature mode. The **W** symbol appears in the display.
- Press and hold the **I** or **Stirn** button and point the sensor tip towards the object or liquid to be measured (never immerse in liquids), maintaining a distance of 3 cm.
- Release the **I** or **Stirn** button. The end of the measuring time is signalled with a short beep and the measured value appears on the display.

Please note that the displayed temperature is the measured and not the adapted surface temperature. It is not the same as the forehead/ear temperature.

To switch back to the Ear/Forehead thermometer mode, press the **M** and **I** buttons simultaneously for 3 seconds until the **W** symbol disappears and there is a short beep. By switching the thermometer off and on again, it automatically exits the Object temperature mode.

The temperature last measured, i.e. only the last value of a series of measurements, is automatically stored as soon as the thermometer switches itself off or is switched off. There are 10 memory spaces available for this.

To switch the thermometer off, press the  and  buttons simultaneously for approx. 3 seconds.

If you do not press any other buttons, the thermometer automatically switches off approx. one minute after displaying the measured value.

To call up the saved values, switch the thermometer on and press the  button.

The date, time, storage space number, the measurement mode symbol and the temperature are displayed in turn.

Browse through the stored measurements by pressing the  button repeatedly.

8 Replacing the battery

When the batteries get weak, the battery warning symbol appears . It is still possible to measure temperature. The batteries must be replaced. When the battery symbol flashes  and  appears in the display, the batteries must be replaced. If the batteries are too flat, then the thermometer will switch off automatically.

Note:

- When changing the batteries, use batteries of the same type, make and capacity.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Use batteries free from heavy metals.

1. Loosen the screw on the battery compartment lid and pull the lid off downwards.
2. Remove the used batteries and insert two new batteries in the direction indicated.
3. Place the battery compartment lid back on and screw it closed.
Dispose of used batteries in accordance with the applicable legal regulations. Never dispose of batteries in the normal household waste.

9 Storage and maintenance

- Clean the sensor tip after each use. Use a clean cloth or cotton bud that can be moistened with disinfectant, alcohol or warm water.
- To clean the entire device, please use a soft cloth slightly moistened with a mild soapy solution.
- Do not use any harsh cleaning products.
- Always store the thermometer with the forehead/cover cap on.
- If prolonged storage is intended, you should remove the battery.
- The device must not be stored or used at an excessively high or low temperature or humidity (see technical data), in sunlight, in association with an electrical current or in dusty locations. Otherwise inaccuracies can occur.

10 Disposal


Batteries should not be placed in the normal household waste.

As a consumer, you are required by law to recycle used batteries. You can recycle your old batteries at public collection points in your community or wherever the same batteries are sold.

Observe the local regulations for material disposal.

Note:

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

Dispose of the unit in accordance with EC Directive **2002/96/EC – WEEE** (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

11 Technical data

Tip: If the device is used other than in accordance with the specifications, perfect functioning cannot be guaranteed. We reserve the right to make technical changes to improve and develop the product.

The accuracy of this thermometer has been carefully checked and developed with regard to a long useful life.

If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means.

Precise instructions for checking accuracy may be requested from the service address.

Name and model	SFT65
Measurement range	Ear/Forehead thermometer mode: 93.2 °F – 109.4 °F (34 °C – 43 °C) Object temperature mode: 32 °F – 212 °F (0 °C – 100 °C)
Laboratory measurement precision	Ear thermometer mode: ±0.4 °F (±0.2 °C) from 95.9 °F – 107.6 °F (35.5 °C – 42 °C), outside of this measurement range ±0.5 °F (±0.3 °C) Forehead thermometer mode: ±0.5 °F (±0.3 °C) from 93.2 °F – 109.4 °F (34 °C – 43 °C) Object temperature mode: ±2.7 °F (±1.5 °C) at < 86 °F (<30 °C); ±5% at ≥ 86 °F (≥ 30 °C)
Interval between 2 measurements	At least 5 seconds
Clinical repeat precision	Ear: Children, 1–5 years: ±0.14 °F (±0.08 °C) Adults: ±0.13 °F (±0.07 °C) Forehead: Children, 1–5 years: ±0.13 °F (±0.07 °C) Adults: ±0.14 °F (±0.08 °C)
Measurement units	°Celsius (°C) or °Fahrenheit (°F)
Ambient operating conditions	60.8 °F – 95 °F (16 °C – 35 °C) with a relative humidity of up to 85% (non-condensing)
Storage conditions	-13 °F to 131 °F (-25 °C to 55 °C) with a relative humidity of up to 85% (non-condensing)
Dimensions	38.2 x 138 x 46.5 mm
Weight	90 g including battery
Battery	2 x 1,5V AAA (LR03)
Memory	For 10 measurements

Symbol key	Device classification type BF 	Keep dry 
Follow instructions for use 		

12 Trouble-shooting

Display	Cause	Solution
<i>H i</i>	The temperature measured is higher than 1) Ear/Forehead thermometer mode: 109.4 °F (43 °C) 2) Object temperature mode: 212 °F (100 °C)	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your retailer or Customer Services.
<i>Lo</i>	The temperature measured is lower than 1) Ear/Forehead thermometer mode: 93.2 °F (34 °C) 2) Object temperature mode: 32 °F (0 °C)	Operate the thermometer only between the specified temperature ranges.
<i>Err</i>	The operating temperature is not in the range from 60.8 °F – 95 °F (16 °C – 35 °C)	Operate the thermometer only between the specified temperature ranges.

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ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1

Guidance and manufacturer's declaration-electromagnetic emissions

The SFT65 Multi-function Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the SFT65 Multi-function Thermometer should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The SFT65 Multi-function Thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The SFT65 Multi-function Thermometer is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Table 2*Guidance and manufacturer's declaration – electromagnetic immunity*

The SFT65 Multi-function Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the SFT65 Multi-function Thermometer should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV air for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T for 0.5 cycle) 40% U_T (60% dip in U_T for 5 cycles) 70% U_T (30% dip in U_T for 25 cycles) 5% U_T (>95% dip in U_T for 5 sec)	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the use of the SFT65 Multi-function Thermometer requires continued operation during power mains interruptions, it is recommended that the SFT65 Multi-function Thermometer be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Table 3*Guidance and manufacturer's declaration – electromagnetic immunity*

The SFT65 Multi-function Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the SFT65 Multi-function Thermometer should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3Vrms 150 kHz to 80 MHz	3Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the SFT65 Multi-function Thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.5 GHz	3V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations from radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can not be predicted theoretic call with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SFT65 Multi-function Thermometer is used exceeds the applicable RF compliance level above, the SFT65 Multi-function Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SFT65 Multi-function Thermometer.

Over the frequency range 150kHz to 80MHz, field strengths should be less than [V1] V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM-for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the SFT65 Multi-function Thermometer

The SFT65 Multi-function Thermometer is intended for use in an electromagnetic environment in which radiated RF distances are controlled. The customer or the user of the T SFT65 Multi-function Thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SFT65 Multi-function Thermometer as recommended below according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 8- MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.