# Spengler

## CARDIOMATE PC<sup>®</sup> & CARDIOMATE PC<sup>®</sup> BLUETOOTH

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- ID: [Text Input]
- Admission No.: [Text Input]
- First Name: [Text Input]
- Outpatient No.: [Text Input]
- Last Name: [Text Input]
- Department: [Dropdown Menu]
- Age: [Text Input]
- Height/Weight: [Text Input] / [Text Input] (with units cm and kg)
- Birthday: [Date Picker] (showing 2020年 6月 5日)
- Systolic/Diastolic: [Text Input] / [Text Input] (with units mmHg)
- Gender: [Dropdown Menu] (showing Others)
- Checked by: [Text Input]
- Audited by: [Text Input]
- Buttons: Previous Patient, Query, OK, Cancel, Create GDT6300, Read GDT6301, Read GDT6302

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FR

**Société/Nom** : Wellkang Ltd

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Enterprise Hub, NW Business Complex, 1 Beraghmore Rd., Derry, BT48 8SE,  
Northern Ireland, UK

**Tél.** : +44(20)3287 6300, +44(33)3303 1126

**Fax** : +44(20)76811874

**Web** : [www.wellkang.ltd.uk](http://www.wellkang.ltd.uk) ; [www.CE-marking.eu](http://www.CE-marking.eu)

### Statement

All contents of this manual are protected by the copyright law. Without the written consent of Shenzhen ECGMAC Medical Electronics Co., Ltd. (hereinafter referred to as «ECGMAC»), any materials contained in this manual shall not be printed, photocopied, reproduced or translated into other languages. ECGMAC owns the copyright and the final interpretation right of this manual.

This manual contains proprietary materials protected by copyright law, including but not limited to technical secrets, patent information and other trade secrets. Users have the obligation of confidentiality and shall not disclose any content in this manual to any third party without written consent. The user's holding of this product manual does not mean that ECGMAC is willing to license the intellectual property rights contained in the product.

The pictures in this manual are interface diagrams and panels for reference only. Please refer to the object. ECGMAC holds the right to revise the contents of each language version of this manual without prior notice.

The product shall be used strictly complying with this manual. User's operation failing to comply with this manual may result in malfunction or accident for which Shenzhen ECGMAC Medical Electronics Co., Ltd. (hereinafter referred to as «ECGMAC») can not be held liable.

### Responsibility of the Manufacturer

When using the product, please ensure that the corresponding electrical installation environment complies with the national standards, and please follow the operation instructions of this manual. ECGMAC only considers itself responsible for any effect on safety, reliability and performance of the equipment if:

- Please use the product according to the product manual;
- The installation, maintenance and upgrading of products shall be carried out by personnel approved or authorized by the ECGMAC;
- The storage environment, working environment and electrical environment of the product should meet the product specifications;
- The label or manufacturing standard of the product is clear and distinguishable, which can prove that the product was manufactured by ECGMAC;
- Machine damage caused by non-human factors such as typhoon, earthquake and other force majeure factors.

Imported and distributed by SPENGLER SAS  
30 rue Jean de Guirmand - 13290 Aix en Provence - France

#### Shenzhen ECGMAC Medical Electronics Co., Ltd.

2nd Floor of Block 2, Haoye Industrial Park, Tiegang Road, Xixiang Street, Baoan District,  
518102 Shenzhen, China

#### Well Kang Ltd

Entreprise Hub, NW Business Complex, 1 Beraghmore Rd, Derry, BT54 8SE, Northern Ireland, UK

<b>1. Safety Guidance</b> .....	<b>68</b>
1.1 Safety Warning .....	68
1.2 Precautions for Use .....	69
1.3 Safety Warnings for Battery Usage.....	70
1.4 List of Symbols .....	70
<b>2.Product Introduction</b> .....	<b>71</b>
2.1 Introduction .....	71
2.2 Product Structure.....	71
2.3 Product Partial Panel .....	73
2.4 Definition for Plug Pins .....	74
2.5 Patient Cable and Electrodes .....	75
<b>3. Before Operation</b> .....	<b>76</b>
3.1 Preparing the Patient.....	76
3.2 Patient Cable Connection .....	76
3.3 Connection of Electrodes .....	77
3.4 Connection between PC and ECG Acquisition Box.....	78
<b>4. Instructions for ECG Acquisition and Analysis Software for Windows PC</b> .....	<b>81</b>
4.1 ECG Acquisition and Analysis Software for Windows PC Installation.....	81
4.2 Main Interface Function Introduction.....	82
4.3 ECG Acquisition and Report Print .....	86
4.4 Reanalysis .....	88
4.5 Patient Record Management.....	89
4.6 System Settings.....	91
<b>5. Instructions for ECG Acquisition and Analysis Software for Android Tablet</b> .....	<b>94</b>
5.1 ECG Acquisition and Analysis Software for Android Tablet Installation.....	94
5.2 Main Interface Function Introduction.....	95
5.3 ECG Acquisition and Report Print .....	96
5.4 Database .....	100
5.5 System Settings .....	100
<b>6. Cleaning, Disinfection and Maintenance</b> .....	<b>106</b>
6.1 General .....	106
6.2 Cleaning.....	106
6.3 Disinfection.....	107
6.4 Care and Maintenance .....	107
<b>7. Common Troubleshooting</b> .....	<b>108</b>
7.1 Common Troubleshooting of ECG Acquisition and Analysis Software for Windows PC.....	108
7.2 Common Troubleshooting of ECG Acquisition and Analysis Software for Android Tablet.....	111
<b>8. Warranty and After-sale Service</b> .....	<b>113</b>
8.1 Warranty .....	113
8.2 Responsibility of Exemption .....	113
8.3 After-sale Service.....	113
8.4 Service Life and Production Date.....	113
<b>Appendix A: Package and Accessories</b> .....	<b>114</b>
<b>Appendix B: Technical Specifications</b> .....	<b>115</b>
<b>Appendix C: Key Components</b> .....	<b>119</b>
<b>Appendix D: List of Functional Accessories</b> .....	<b>119</b>
<b>Appendix E: Network Security Specifications</b> .....	<b>119</b>
<b>Appendix F: EMC Information</b> .....	<b>121</b>
Table 1.....	121
Table 2.....	122
Table 3.....	122
<b>Appendix G: Manufacturer Information</b> .....	<b>123</b>
<b>Appendix H: EU Representative Information</b> .....	<b>123</b>

# **E 1. Safety Guidance**

## **1.1 Safety Warning**

- ⚠ This instrument is suitable for the detection and recording of routine human ECG.
- ⚠ This system is not designed for intracardiac use or direct cardiac application.
- ⚠ This instrument is not intended for home use.
- ⚠ This instrument is not intended for treatment or monitoring.
- ⚠ This instrument is provided for the use of qualified physicians or personnel professionally trained.
- ⚠ The operator is supposed to be familiar with the contents of this User Manual before operation.
- ⚠ Do not use this instrument near anesthetic gases, oxygen, hydrogen or other flammable and corrosive chemicals, otherwise there will be a risk of explosion.
- ⚠ Do not use this instrument in the hyperbaric oxygen chamber, otherwise there will be a risk of explosion.
- ⚠ Do not use this instrument near high-voltage, high- static, otherwise it may discharge instantly and produce sparks.
- ⚠ Only service engineers authorized by the manufacturer can open the shell.
- ⚠ Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC/EN standards (e.g. IEC/EN 950 for data processing equipment and EN 60601-1 for medical equipment). Furthermore all configurations shall comply with the valid version of the standard EN 60601-1-1. Therefore anybody, who connects additional equipment to the signal input or output connector to configure a medical system, must make sure that it complies with the requirements of the valid version of the system standard EN 60601-1-1. If in doubt, please consult our service department or the local agent on your side.
- ⚠ When defibrillator is used simultaneously with the instrument, the operator should not touch patient and bed. All the electrodes (whether connected to the patient or not) and the patient do not need to be grounded. AgCl coated chest electrode and patient cable with the function of defibrillation should be selected. If the defibrillation time is more than 5 seconds, disposable chest electrode should be used to avoid skin burn of patients by metal electrode. Please do not use it at the same time with other electrical stimulation equipment. If it is necessary to use them at the same time, professional and technical personnel shall be present for guidance, and the selected accessories shall be the accessories designated by the manufacturer.
- ⚠ Be careful when patient is connected to multiple equipments, the total leakage current might cause the injury. Only equipment compliant with IEC60601-1:2006 is allowed to connect with this instrument, when it is connected with other equipments. Meanwhile, reliable connection between the potential equalization needs to be carefully considered. After the connection, the users must measure the total leakage current by themselves to determine whether it meets the requirements or usage condition.
- ⚠ The pacemaker (if installed) may influence the accuracy and analysis result. Under this situation, it is suggested that the doctors should identify and analyze it according to the waveform. In addition, it will also increase the potential risk. In this case, when recording ECG, pay special attention to safety issues, and appropriate measures must be taken to ensure that the leakage current is at a safe level.
- ⚠ This instrument can not be used with high-frequency electrotome.
- ⚠ The operator should not touch patient, patient bed, table or instrument when it is used simultaneously with the defibrillator or pacemaker.
- ⚠ Please use the patient cable and other accessories supplied by ECGMAC. Otherwise, it will affect the performance and function, even damage the ECG instrument.
- ⚠ Be sure that all electrodes have been connected to the correct position of human body before operation. Avoid electrodes (including neutral electrodes) and patient to be contacted with any conducting parts or earth.
- ⚠ This instrument is equipped with corresponding user permission. Please properly keep the authorized account and password information, and do not disclose the account and password information to others.
- ⚠ The patient's health information can only be used by authorized medical personnel for cardiac function diagnosis and analysis, and shall not be disclosed to a third party.


## 1.2 Precautions for Use

- ❗ The instrument should be placed in a horizontal table to avoid excessive vibration and shock during movement.
- ❗ Avoid liquid splash and excessive temperature. The temperature must be kept between +5°C to + 40°C while working.
- ❗ Do not use the equipment in dusty environment with poor ventilation or in the presence of corrosive, such as salt, sulphur and chemicals.
- ❗ Be sure that there is no intense electromagnetic interference source around the equipment, such as radio transmitter or mobile phone, electrosurgical equipment, ultrasound equipment, radiological equipment and magnetic resonance imaging equipment etc.
- ❗ Before using the instrument, it is necessary to check whether there is any damage to the instrument, patient cable and electrodes that may affect the safety of the patient. Replacement should be taken if there is any evident defectiveness or aging symptom which may impair the safety or performance.
- ❗ In order to record ECG data accurately, the instrument should be placed in a quiet and comfortable environment.
- ❗ The equipment and reusable accessories can be sent back to the manufacturer for recycling or proper disposal after their useful lives.
- ❗ If any accident happens during use, please immediately stop using the machine.
- ❗ Please make sure that there are no other Bluetooth devices open in the surrounding environment except for the computer installed with ECG acquisition system software and ECG acquisition box when transferring data via Bluetooth of the Bluetooth ECG.
- ❗ It is forbidden to use the mobile phone in the surrounding environment to avoid data error when transferring data via Bluetooth of the Bluetooth ECG.
- ❗ Do not touch or move the acquisition box or USB connection cable when transferring data via USB cable. It is recommended that the USB cable should not exceed 1.5m, otherwise data loss may occur.
- ❗ Do not move the acquisition box when transferring data via Bluetooth. The barrier free transmission range between ECG acquisition and tablet or PC shall not exceed 5 meters, otherwise data loss may occur.
- ❗ In case of data pause during Bluetooth transmission, please check the connection distance between the acquisition box and the software or the surrounding environment.
- ❗ Please choose to overlay installation when upgrading the software, otherwise the data will be lost.
- ❗ The software of ECG signal acquisition system can only receive the import of ECG data with correct format.
- ❗ In the process of using, the instrument will save the patient's patient record information and ECG data records. It is recommended to regularly backup the data to prevent the loss of data caused by computer damage.
- ❗ Please close other applications before using this instrument. Please do not use other software on the same computer while running this program.
- ❗ It is suggested that a specific computer should be used with this instrument, and irrelevant software should not be installed.
- ❗ The USB flash disk / SD card used to import and export ECG data needs to be scanned by antivirus software.
- ❗ When using the instrument, the computer screen protection, automatic sleep and sleep functions shall be disabled.
- ❗ This instrument only provides ECG as the reference for medical staff. The data itself does not have the significance of diagnosis results, and the final diagnosis should be provided by medical staff.

## 1.3 Safety Warnings for Battery Usage

- ⚠ The battery is to be disposed of according to local regulations after their useful lives.
- ⚠ Do not ruin the battery by hammering, beating, chiseling a metal into it or other else, which will cause the battery to be deformation, heat, to smoke, burn and other danger.
- ⚠ When leakage or foul smell found, stop using the battery immediately. If the leakage liquid gets to your skin or cloth, cleanse it with clean water at once. If the leakage liquid gets into your eyes, do not wipe them. Irrigate them with clean water first and go to see a doctor immediately.
- ⚠ When the battery expires, or foul smell, leakage, deformation, discoloration and distortion have been found, stop using it.
- ⚠ Take out the battery when not using the equipment for long time.
- ⚠ Do not reverse the anode and the cathode when installing the battery, otherwise it may cause explosion.
- ⚠ Children are not allowed to touch the battery.

## 1.4 List of Symbols

Symbols	Function	Symbols	Function
	USB Port		Bluetooth
	Lot Number		Product Serial Number
	Manufacturer		Date of Manufacture
	Recycle		Transportation and Storage Temperature Range
	Avoid Direct Light		Keep Dry
	Tier Limitation		DEFIBRILLATION-PROOF TYPE CF APPLIED PAR
	Fragile		Note (general warning): The Information You Should Know to Avoid Possible Damage to Patients or Operators.
	Note: The Information You Should Know to Avoid Possible Damage to The Equipment		Transportation and Storage Atmospheric Pressure
	Transportation and Storage Humidity Range		Follow Instructions for Use
	CE		European Authorized Representative
	Recycling electronic equipment		

**Note :** The «acquisition box» in this manual is a general designation of Bluetooth ECG acquisition box and PC-ECG acquisition box.

## 2. Product Introduction

### 2.1 Introduction

PE-12 series PC-ECG and Bluetooth ECG (hereinafter referred to as PC-ECG and Bluetooth ECG) acquire ECG signals from adults or children under resting status through body surface electrodes. The recorded ECG waveform and the measurement parameters and interpretation provided by the instrument serve as an important basis for the diagnosis of heart disease. PE-12 series products are light and compact. They support 12 leads and 6 leads acquisition methods, multiple display and report formats, PDF, JPEG, BMP, EM-XML, FDA-XML, SCP, DICOM, BKG and other formats for data report, and ECG waveform replay. PE-1201 series products are PC-ECG including PE-1201 and PE-1201A. PE-1204 serials products are Bluetooth ECG including PE-1204 and PE-1204A. PE-1204 and PE-1204A support USB and Bluetooth, and can support Windows operating system and Android operating system. PE-1201 and PE-1201A only support USB cable acquisition mode, and can support Windows operating system or Android operating system.

**Scope of application:** it is used for medical institutions to extract ECG wave of human and generate ECG measurement parameters and interpretation to diagnose heart diseases and facilitate the relative researches.

The main use occasions include hospitals or medical institutions. The PC-ECG and Bluetooth ECG should be used by doctors or trained medical professionals in hospitals or medical institutions.


**Contraindications:** None

### 2.2 Product Structure


The PC-ECG and Bluetooth ECG system include ECG Acquisition Box, Patient Cable, USB cable, Dongle, Chest Electrodes, Limb Electrodes and ECG acquisition and analysis software.

#### 2.2.1 PC-ECG/Bluetooth ECG Connection Sketch-map

##### 2.2.1.1 Sketch-map of Connecting Computer and Acquisition Box via USB

##### 2.2.1.2 Sketch-map of Connecting Bluetooth ECG Acquisition Box and Computer (or Android Tablet) via Bluetooth (only applicable to PE-1204, PE-1204A)

## NE 2.2.2 Product Model and Specification

Model	Voltage of Power Supply	Sampling Rate	Number of Leads	Configuration	Dimension (mm)	Difference
PE-1201	d.c.5V	2000	12 leads	Main unit (ECG acquisition box) ECG acquisition and analysis software Dongle (USB Dongle or License Code) ECG patient cable USB cable Chest electrodes and limb electrodes (optional)	126×68×24	No Bluetooth module Battery power not supported Casing Color: Warm Gray 12 bit ADC 10-lead Patient Cable
PE-1204	d.c.5V ou d.c.3V	2000	12 leads	Main unit (Bluetooth ECG acquisition box) ECG acquisition and analysis software Dongle (USB Dongle or License Code) Patient cable USB cable Chest electrodes and limb electrodes (optional)	126×68×24	4.0 Bluetooth module Power supplied by AAA alkaline battery Casing Color: Warm Gray 12 bit ADC 10-lead Patient Cable
PE-1201A	d.c.5V	32000	12 leads	Main unit (ECG acquisition box) ECG acquisition and analysis software Dongle (USB Dongle or License Code) ECG patient cable USB cable Chest electrodes and limb electrodes (optional)	126×68×24	No Bluetooth module Battery power not supported Casing Color: Warm Gray 24 bit ADC 10-lead Patient Cable
PE-1204A	d.c.5V ou d.c.3V	32000	12 leads	Main unit (Bluetooth ECG acquisition box) ECG acquisition and analysis software Dongle (USB Dongle or License Code) Patient cable USB cable Chest electrodes and limb electrodes (Optional)	126×68×24	4.0 Bluetooth module Power supplied by AAA alkaline battery Casing Color: Warm Gray 24 bit ADC 10-lead Patient Cable

## 2.3 Product Partial Panel


### 2.3.1 Main Unit Front Panel and Back Panel

Sign	Name	Description
A	USB Port	Used for connecting USB data cable to USB port of computer
B	Patient Cable Socket	Used for collecting ECG patient cable, the other end of it is connected with the electrodes


### 2.3.2 Main Unit Top View

#### 2.3.2.1 PE-1204/PE-1204A Main Unit Top View

Sign	Name	Description
C	Power	When the device is turned on, the light is on
D	Bluetooth	This light flashes when the device transmits data in Bluetooth mode.
E	USB	This light flashes when the device transmits data via USB.


#### 2.3.2.2 PE-1201/PE-1201A Main Unit Top View

Sign	Name	Description
G	Power	When the device is turned on, the light is on
H	Work	This light flashes when the device transmits data via USB.
K	Electrode Check	If the lead falls off, the indicator light will be on.


**2.3.3 Main Unit Bottom Panel**

*2.3.3.1 PE-1204/PE-1204A Main Unit Bottom Panel*

Sign	Name	Description
F	Battery Compartment	Battery compartment cover (please install 2 AAA alkaline batteries)


*2.3.3.2 PE-1201/PE-1201A Main Unit Bottom Panel*

Sign	Name	Description
F	Empty compartment	No function definition

**2.4 Definition for Plug Pins**

**2.4.1 Patient Cable Socket**


 Applied part of type CF with defibrillator proof

### 2.4.2 Definition for Plug Pins

Pin	Signal	Pin	Signal	Pin	Signal
1	C2 (input)	6	SH	11	F (input)
2	C3 (input)	7	NC	12	NC
3	C4 (input)	8	NC	13	C1 (input)
4	C5 (input)	9	R (input)	14	NC
5	C6 (input)	10	L (input)	15	N or RF (input)


## 2.5 Patient Cable and Electrodes

### 2.5.1 Patient Cable


Patient cable includes main cable and lead wires. The lead wires include 6 chest lead wires and 4 limb lead wires. The user can identify them through the colors and the identifiers on the connectors.

### 2.5.2 Chest Electrodes

## EN 2.5.3 Limb Electrodes


### 3. Before Operation

#### 3.1 Preparing the Patient

1. For the first time patient, the operator shall communicate and explain in advance:
  - 1) Do not smoke or drink tea, coffee or alcohol for 2 hours before examination;
  - 2) The patient should wear loose clothes as much as possible to facilitate ECG examination;
  - 3) Before receiving ECG examination, the patient should rest for a while, keep calm and avoid being nervous;
  - 4) If there is dirt on the skin where the electrode is placed, the skin should be cleaned first;
  - 5) If there is too much hair where the electrode placed, the hair should be shaved to reduce the resistance.
2. When possible, prepare the patient in a quiet room or area where others can't see the patient. If there are other people in the room, the curtain should be pulled at the bedside during ECG acquisition.
3. Ensure that the patient can lie down comfortably. The more relaxed the patient is, the less interference the ECG wave will suffer.
4. Once the electrode and patient cable are connected, inform the patient: breathe naturally, keep still, don't talk, don't chew or clench the teeth.

#### 3.2 Patient Cable Connection

Plug the patient cable connector into Patient Cable Socket on acquisition box as shown below, and lock the screw on both sides.

- ❗ Only use the original patient cable of manufacturer equipped. Do not use any other models.
- ❗ Patient cables socket is designed for ECG signal input only. Do not use it for other Purpose.

### 3.3 Connection of Electrodes

The installation of electrodes is critical for accurate record of ECG signals. Please ensure electrodes are well-connected. Disposable electrodes cannot be used together with reusable electrodes (old or new). Different models cannot be used together otherwise it will seriously affect the accuracy of ECG record. Electrodes or lead wire socket cannot touch other conductors, such as metal beds. All electrodes need to be replaced together during replacement.

#### 3.3.1 Connection of Limb Electrodes

Limb electrodes should be installed on the soft skin of hands and legs. Clean the skin where electrodes will be installed with medical alcohol; Apply a small amount of electrode gel on the skin. The connection is shown in the right figure.


American Label	European Label	Electrode Placement
RA	R	Right arm
LA	L	Left arm
RL	N or RF	Right leg
LL	F	Left leg

#### 3.3.2 Connection of Chest Electrodes

Clean the skin where electrodes will be connected with alcohol (75%). Apply a small amount of electrode gel to skin where had been cleaned and the edge of the suction ball of the chest electrodes. Press the suction balls of the chest electrodes and attach them to C1 to C6. The connection is shown in the following figure:

- C1 (V1): Fourth inter-costal space at right border of the sternum.
- C2 (V2): Fourth inter-costal space at left border of the sternum.
- C3 (V3): Midpoint between C2 and C4.
- C4 (V4): The intersection of the fifth inter-costal space of the left sternum and the central line of the clavicle.
- C5 (V5): Left anterior axillary line on the same horizontal level as V4.
- C6 (V6): Left mid-axillary line on the same horizontal level as V4.

#### ⚠ Note:

Electrode gel coating should to be separated from each other. Chest electrodes should not touch with each other to avoid short circuit. If no electrode gel is available, apply a small amount of alcohol (75%) instead. Clean the electrode placement site with alcohol (75%), and then install the corresponding electrode to ensure that the skin of electrode contact area is moderately moist. Do not use physiological saline to replace electrode gel. Otherwise, it will cause electrode corrosion.

### 3.3.3 Electrode Identification and Connection Color Table

Electrode Position	IEC		AHA	
	Electrode identifier	Electrode color code	Electrode identifier	Electrode color code
Right arm	R	Red	RA	White
Left arm	L	Yellow	LA	Black
Right leg	N or RF	Black	RL	Green
Left leg	F	Green	LL	Red
Chest 1	C1	Red	V1	Red
Chest2	C2	Yellow	V2	Yellow
Chest3	C3	Green	V3	Green
Chest4	C4	Brown	V4	Blue
Chest5	C5	Black	V5	Orange
Chest6	C6	Purple	V6	Purple

## 3.4 Connection between PC and ECG Acquisition Box


### 3.4.1 Connect the PC and the Acquisition Box via USB Cable

#### 3.4.1.1 USB Cable Connecting PC and Acquisition Box

1) Connect the attached USB data cable to the computer and the acquisition box (as shown below).

2) « System Settings » → « Device and Sample settings » → « Device » of ECG acquisition and analysis software for Windows PC is set as «USB ECG Device».

### 3.4.1.2 Connect Android Tablet and Acquisition Box via USB Cable

Android Tablet (hereinafter referred to as tablet) can only be connected with ECG acquisition box via USB cable and OTG cable at the same time. The USB Master Port is connected with the acquisition box, the standard port of the USB cable is connected with the USB wide Port of the OTG cable, and the MINI USB port of the OTG cable is connected with the tablet. After the connection between the acquisition box and the tablet is completed, open the ECG acquisition and analysis software for Android Tablet, click **【Sampling】** for the first time, and the following dialog box will pop up:

Select «Use by default for this USB device» and select **【OK】** to obtain the USB access authority successfully.

**! Note:** please choose the OTG cable model recommended by the manufacturer.

### 3.4.2 Connect PC and Bluetooth ECG Acquisition Box through Bluetooth

Bluetooth wireless acquisition mode is only applicable to PE-1204 and PE-1204A.

#### 3.4.2.1 Use Bluetooth to Connect PC and Bluetooth ECG Acquisition Box

The computer installed with ECG acquisition and analysis software for Windows PC is connected to the Bluetooth ECG acquisition box through Bluetooth, and the steps are as follows:


- 1) Open the battery compartment cover of the Bluetooth ECG acquisition box, correctly install two AAA alkaline batteries, and then close the battery compartment cover.
- 2) Turn on the Bluetooth function of the computer or connect the Bluetooth adapter to the computer (please choose the Bluetooth adapter model recommended by the manufacturer).
- 3) Press and hold the power button of the acquisition box, and release it when the power light, Bluetooth connection status light and USB connection status light are all on.
- 4) Select Bluetooth on the computer, and the Bluetooth name of the Bluetooth ECG acquisition box is «ECGMAC». Correctly connect them and ensure the successful Bluetooth pairing.
- 5) In the «Bluetooth and other settings» interface of the computer, click «more Bluetooth options - COM port», click «COM port», and check the corresponding COM number of ECGMAC.
- 6) Select Bluetooth ECG Device in «Settings» - >»Device and Sample Settings» - >»Device» - > of ECG acquisition and analysis software for Windows PC, set COM NO. The same with the COM number of ECGMAC. The software then can obtain data from wireless signal through Bluetooth transmission protocol.

**! Note :**


- 1) The driver of Bluetooth adapter must be installed correctly.
- 2) When the Bluetooth ECG acquisition box is connected to PC via Bluetooth, you must choose «Device» ->»Bluetooth ECG Device» and accurately input the Bluetooth outgoing COM number randomly assigned by the computer.
- 3) When Bluetooth is connected, the Bluetooth ECG acquisition box only supports power supply from two AAA Alkaline battery. Please check that the batteries are installed correctly.

**EN** 3.4.2.2 Connection between Android Tablet and Bluetooth ECG Acquisition Box via Bluetooth

- 1) Open the battery compartment cover of the Bluetooth ECG acquisition box, correctly install two AAA alkaline batteries, and then close the battery compartment cover. Press and hold the power key of the Bluetooth ECG acquisition box, and release it when the power light, Bluetooth connection status light and USB connection status light are all on.
- 2) Open Tablet ECG, and click **【Sampling】**. The following prompt box will pop up:

- 3) Select **【OK】** to pair with the opened Bluetooth ECG acquisition box. Then, automatically acquisition of ECG waveform after successful pairing will start.

## 4. Instructions for ECG Acquisition and Analysis Software for Windows PC


### 4.1 ECG Acquisition and Analysis Software for Windows PC Installation

#### 4.1.1 System Installation and Operation Environment



Configuration	Specific Requirement
CPU	≥1.2 G
System Memory	512MB or above
Main Board	Inter main board recommended
Hard Disk	40GB or above
Display card	Over 1m
Display	≥43 cm(17 Inch), >1024*768, ≥16
CD-ROM	Above 24 times speed
Operating System	Microsoft Windows 2000/xp/7/8/10
Printer	Resolution above 600dpi
Others	Bluetooth 4.0

#### 4.1.2 Installation du 4.1.2 Software System Installation logiciel


- Put the installation CD into the CD-ROM, and the software will automatically start the installation (if it cannot be installed automatically, please execute the CD root directory SetupECGApp.exe File).

- Click **Next** ; choose program installation path in the “choose destination location”, and click **Next** .

- EN**
- 3) In the «ready to install» prompt box of the installation program, click **【Install】** to start the installation process.
  - 4) The installation program will automatically install the driver of the acquisition box and the driver of the Dongle.
  - 5) After the installation of Dongle driver, the interface for installing USB acquisition box driver will pop up, and click **【Install】** to start the installation of acquisition box driver.
  - 6) After installation, in the «InstallShield Wizard» of the installation program, click **【Finish】**.

**! Note :** The software installation will automatically install the Driver for the Dongle and USB acquisition box. If the driver installation fails, please install the driver manually.


## 4.2 Main Interface Function Introduction

### 4.2.1 System Startup

1. Turn on the the computer and printer.
2. Insert the USB Dongle into the USB port of the computer installed with ECG acquisition and analysis software for Windows PC.
3. Run the software as an administrator and enter the ECG acquisition and analysis software for Windows PC initialization main interface (as shown below):

Method 1: right click the software icon and choose to run as administrator. Follow this procedure every time you open the software.

Method 2: right click the software icon and select Property → Compatibility → Run as administrator.


**Note:** Please make the following checks before starting the ECG acquisition and analysis software for Windows PC


- 1) Check whether the ECG acquisition and analysis software for Windows PC has been installed.
- 2) Check whether the Dongle is connected effectively.

### 4.2.2 Main Interface and Toolbar Introduction

The main interface display includes acquisition interface and ECG waveform replay interface. The acquisition interface displays the acquired ECG waveform in real time. After the acquisition, it enters the replay interface. Through moving the progress bar, the acquired ECG waveform in the whole process can be viewed back.







Sampling


Replay

As shown in the above figure, the toolbar contains different buttons. See the following table for specific functions:

Icon	Name	Function
	Sampling	Click <b>【Sampling】</b> , input the corresponding patient information in the pop-up “patient information” dialog box, and click <b>【OK】</b> to start acquiring ECG data.
	Stop	You can click “Stop” to stop the acquisition.
	Paper speed	Switch the paper speed of the current main interface. Automatically switch paper speed: click <b>【Paper speed】</b> to automatically cycle paper speed. Manually switch paper speed: click the <b>【Paper speed】</b> drop-down box to manually select paper speed, where “√” represents the paper speed of the current interface.
	Gain	Switch the gain of the current main interface Automatically switch gain: click <b>【Gain】</b> to automatically cycle the gain. Manually switch gain: click the drop-down box of <b>【Gain】</b> to manually select the gain, where “√” represents the gain of the current interface. Note: when the gain is set to 10 / 5 mm / mV, 20 / 10 mm / mV and AGC (automatic gain), the reference gain of the main interface is displayed as 10 mm / mV, and the gain of the channel is not 10 mm / mV. The multiple relationship is added beside the lead symbol, where (x2), (x4), (x1 / 2), (x1 / 4) and (x1 / 8) represent the gain of the current channel as 20 mm / mV, 40 mm / mV, 5 mm / mV, 2.5 mm / mV and 1.25 mm / mV respectively.
	Filter	The filter parameter value can be adjusted during the acquisition. Filter includes AC filter, EMG filter, Baseline filter and Low pass filter.
	Upload (optional)	According to the upload method set by the system, upload the current patient record report to the corresponding server.
	Format	Quickly switch the lead display format of the acquisition interface. Automatically switch format: click <b>【Format】</b> to switch the display format automatically. Manually switch format: click the <b>【Format】</b> of drop-down box to select the format manually, where “√” is the display format of the current interface.
	Database	Enter the database interface to query, edit, import, export and upload the patient record.
	Reanalysis	Reanalyze the current replayed ECG data. The parameters and conclusions can be manually modified. Please refer to [4.4 Reanalysis] for details.
	Settings	Click <b>【Settings】</b> to enter the system setting interface, where you can set related information.
	Print	Click <b>【Print】</b> on the toolbar to print the ECG report. Click the drop-down box of <b>【Print】</b> in the toolbar to select the report format.
	Minimize	Click <b>【Minimize】</b> to minimize the software interface, and only keep the icon in the taskbar.
	Exit	Exit the system

### 4.2.3 Right Click Function Introduction

In the replay interface, the right click displays the following functions: export, ruler, parallel ruler, aVR, aVL leads correction, magnifier, print and upload (optional), as shown in the following figure:

Name	Function
Export	The export of data files and image reports. See the export settings in 4.5 for details.
Ruler	Measure the waveform width and amplitude: Press and hold the left mouse button at the starting point where the waveform needs to be measured, drag the mouse to the measuring ending point, release the mouse, and the ruler will display the time difference and amplitude difference between the two positions.
Parallel ruler	You can check the regularity of heart beat rhythm: Right click to select "parallel ruler", click the left mouse button on the waveform of the main interface and drag the mouse to show a series of parallel lines. The software calculates the interval time difference and heart rate value through the distance between the parallel lines.
aVR, aVL leads correction	ECG data of related channels of left and right limb leads can be exchanged: Right click to select "aVR, aVL leads correction" in the main interface, select "yes" to continue aVR, aVL leads correction, and select "no" to give up the operation.
Magnifier	Amplification of ECG waveform: Right click to select "magnifier", press and hold the left mouse button at any position of the waveform, and drag the mouse to enlarge the ECG waveform at the mouse position.
Print	Right click to select "print", and directly print the ECG report in the default print format. If the preview function is on, select Print to display the print preview interface.
Upload (optional)	According to the preset upload setting, upload the currently replayed patient record report to the corresponding server.

## EN 4.3 ECG Acquisition and Report Print

Bluetooth ECG can acquire ECG data through USB cable acquisition and Bluetooth wireless acquisition. Before acquiring ECG data, please make sure that the patient has connected the electrode, the patient cable correctly connected to the ECG acquisition box, and the acquisition box reliably connected to the computer. See [3. Before Operation] for details.

### ⚠ Note :

- 1) Check that the connection mode between the computer and the acquisition box is consistent with the selection of «Settings» → «Device and Sample settings» - >»Device».
- 2) In case of waveform pause or acquisition interruption during acquisition, please check the reliability of connection between acquisition box and software.
- 3) When the waveform interference is large or the lead falls off during the acquisition, please check the reliability of the connection between the patient cable and the acquisition box and the human body.
- 4) In the process of Bluetooth acquisition, if the distance between software and the Bluetooth ECG acquisition box exceeds the effective transmission distance of Bluetooth, the acquisition will be interrupted automatically.

### 4.3.1 Patient Information and ECG Acquisition

1. Open ECG acquisition and analysis software for Windows PC, click **【Sampling】**, and input corresponding information in the pop-up “Patient Information” dialog box.

The screenshot shows a 'Patient Information' dialog box with the following fields and controls:

- ID: [Text Input]
- Admission No.: [Text Input]
- First Name: [Text Input]
- Outpatient No.: [Text Input]
- Last Name: [Text Input]
- Department: [Dropdown Menu]
- Age: [Text Input]
- Height/Weight: [Text Input] / [Text Input] (units: cm / kg)
- Birthday: [Date Picker] (2020年 6月 5日)
- Systolic/Diastolic: [Text Input] / [Text Input] (units: mmHg / mmHg)
- Gender: [Dropdown Menu] (Others)
- Checked by: [Dropdown Menu]
- Audited by: [Dropdown Menu]
- Buttons: Previous Patient, Query, OK, Cancel, Create GDT6300, Read GDT6301, Read GDT6302

The patient information interface includes the following:

Name	Function
Patient Information	You can input patient ID, name, age, birthday, gender, examining doctor, examining doctor, inpatient number, outpatient number, department and user-defined items.
Previous Patient	Click <b>【Previous Patient】</b> to automatically fill in the previous patient information.
Worklist Update (Optional)	<ol style="list-style-type: none"> <li>1. Start worklist query in system settings and download patient information from the specified server.</li> <li>2. Auto: select auto. According to the automatic update time set by the system, automatically query the patient information of the day from worklist server at regular intervals. The update interval can be set.</li> <li>3. Manual: it is required to manually query patient information from worklist server.</li> <li>4. Update range: query the server patient information according to the selected update time range. The update scope includes the range of today / tomorrow / yesterday / date.</li> </ol>
Query	<ol style="list-style-type: none"> <li>1. Enable worklist (Optional) query in the system setting, input the corresponding patient information conditions, and click <b>【Query】</b> to select the patient information data queried from worklist server.</li> <li>2. Enable HL7 (Optional) query in the system settings, input ID, click query to query the patient information of corresponding ID from the server.</li> <li>3. If Worklist and HL7 are disabled, the patient information data saved in the local database is queried according to the input patient information conditions, which is convenient for direct use.</li> </ol>

Name	Function
OK	Save the current input patient information and enter the main interface to start sampling.
Cancel	Do not save the current input information and return to the main interface.

## 2. Start Sampling

In the patient information interface, click **【OK】** to start sampling ECG data.

Presampling: the «Presampling Time» of the “System Settings”- >» Device and Sample Settings «is set to the non close state. ECG data under presampling mode will not be recorded.

After the set presampling time, the system starts to record ECG data, and the recording time of the main interface starts to count.

In the acquisition process, click **【Sampling】** to clear the recording time and restart recording ECG data for the current patient.

It is strongly recommended that when there is strong interference in the current acquisition process, click **【Sampling】** again to eliminate the interference after the interference disappears, so as to obtain high-quality ECG waveform and reduce misdiagnosis.

## 3. Stop

Manual stop: click **【Stop】** to end the current acquisition. When the recording time is less than 10s, the software does not save the data.

Auto stop: if Auto is selected in “Settings” - > “Device and Sample Settings”, “Sample Mode”, ECG acquisition and analysis software for Windows PC will automatically stop and save the file after the record reaches the set sampling time.

If the waveform is replayed after the acquisition is stopped, the quality of the waveform is poor. It is strongly recommended to re-sample the patient to obtain high-quality ECG waveform and reduce misdiagnosis.

### 4.3.2 Check the Quality of ECG

#### 4.3.2.1 Lead off Alarm

During the acquisition process, if the lead connection is unstable or dropped, the software will prompt the alarm of lead off. Check the loose electrode or lead according to the prompt information, and reconnect the electrode or lead.

#### **!** Note:

- 1) When there is a lead off alarm during the acquisition process, it is suggested to acquire again.
- 2) The data acquisition under the condition of lead off may lead to information loss, missed diagnosis or misdiagnosis.

#### 4.3.2.2 Eliminating Waveform Interference

If there is interference in the waveform, please check the filter settings so that the filter settings are most suitable for the current environment.

If the filter cannot eliminate the interference, please check the electrode connection and skin condition of the patient, and check whether the ambient temperature is appropriate.

If the above steps still can not eliminate interference, please check whether there is strong interference equipment in the surrounding environment.

### 4.3.3 Report print

#### 4.3.3.1 Manually print

Manual printing can be realized in the following three ways:

1. Click **【Print】** on the toolbar to print the ECG report in the default print format.

- EN**
2. Select the printing format through the **【Print】** drop-down triangle, and then print the ECG report according to the selected printing format.
  3. Right click the replay interface and select **【Print】** to print the ECG report in the default print format.

#### 4.3.3.2 Auto Print

«System settings» → «print settings» to select automatically print the report. ECG acquisition and analysis software for Windows PC will automatically print the ECG report according to the set conditions after the acquisition stops.


#### 4.3.3.3 Print Preview

The current ECG report can be previewed when printing the report. «System settings» → «Print setting» to choose print preview function when printing». When printing is selected, enter the print preview interface first.

## 4.4 Reanalysis

In the replay interface of ECG waveform, click **【Reanalysis】** to enter the interface of reanalysis, as shown in the figure below. The interface of reanalysis shows ECG characteristics, average template of each lead, interpretation, etc., and ECG parameters can be re-analyzed.

Modify parameters: users can manually input the parameters or move the location line of the average template to modify the parameters.


Switch display channel: users can click the lead name to switch to display the average template of the corresponding lead.

Modify conclusion:

- Add: directly input the conclusion in the conclusion input box manually or select the conclusion in the analysis result index, and click [Add] to add the conclusion.


- Clear: clear the conclusion information.
- Save as template: save the current conclusion as user-defined template.
- Edit template:
  - ◆ Import template: import user-defined templates in bulk.
  - ◆ Update: update the selected template content.
  - ◆ Save as template: save the edited template as a new template.
  - ◆ Delete: delete the selected template.
  - ◆ Close: close the edit template box.

- Reset: restore ECG characteristic parameters and conclusions to the initial state.
- OK: save the modified information and close the reanalysis interface.
- Cancel: do not save the modified information and close the reanalysis interface.

## 4.5 Patient Record Management

In the main interface, click **【Database】** to enter the patient record interface, where you can query, open, delete, edit, compare, upload (optional), import and export patient records. The patient record interface is mainly composed of patient record query, list of patient record and management tools.

### 4.5.1 Patient Record Query

Patient Record Query can quickly locate relevant data by setting relevant conditions. After

**EN** setting the query conditions, click **【Query】** to realize the query function, and count the number of data queried. To cancel the query result display, click **【Clear all conditions】** .

Query conditions	To query according to patient information such as as ID, name, age, gender, birthday, department, check doctor, audited doctor, admission number, outpatient number and user-defined information.
Check date	The query results only display the patient records within the set check time range.
Upload condition	There are two upload results: not uploaded and uploaded. The upload status query only takes effect when "√" is selected.
Leads	It includes three options: all leads, 12 leads and 6 leads. All leads display the data of 12 leads and 6 leads at the same time; 12 leads only display the data of 12 channels; 6 leads only display the data of 6 channels.

#### 4.5.2 List of Patient Record

The patient record list displays all eligible patient records, and displays the check time, ID, name, age, gender, uploaded, number of leads, file path, view status, etc.

Viewed: indicates whether the file is viewed. The default viewing status of the data collected by the software is yes, and the default imported data is No.

Left click	Select a single patient record
Shift + left click	Continuously select multiple patient records
CTRL + left click	Not continuously select multiple patient records
Right click to "Check All"	Select all patient records
Right click to "Uncheck All"	Cancel select all patient records
↑	Click "↑" on the keyboard to select the previous patient record
↓	Click "↓" on the keyboard to select the next patient record
←	To press and hold the "←" key on the keyboard,you can drag the scroll bar of patient record to the right
→	To press and hold the "→" key on the keyboard,you can drag the scroll bar of patient record to the left


#### 4.5.3 Management Tools

Name	Function
Open	Open the selected patient record, enter the main interface, and display the ECG waveform and patient information of the selected patient record.
Delete	Delete the selected patient record.
Edit	Select the patient record needed to be edited and click <b>【Edit】</b> to edit the patient information.
Sampling	Click <b>【Patient record】</b> and click <b>【Sampling】</b> to automatically use the patient information of the patient record for ECG data acquisition again.
Compare	Only when the status of two patient records is selected, the <b>【Compare】</b> function takes effect, and the information of the two selected patient records will be displayed on the same screen. The gain, paper speed, format, measurement and printing functions can also be set in the comparison interface.

Name	Function
Import	Import ECG data from disk, support EM-XML, FDA-XML, SCP, BKG, DICOM standard format ECG data.
Export	Support the export of ECG data from software to disk.
Upload (optional)	Select one or more patient records and click <b>【Upload】</b> to upload the selected patient record data to the corresponding server according to the set upload method. Upload status: <b>Yes</b> indicates uploaded <b>No</b> indicates not uploaded.
Close	Click <b>【Close】</b> or click the close icon «X» in the upper right corner to exit the patient record management interface and return to the main collection interface.

## 4.6 System Settings

Click **【Settings】** to open the administrator password input dialog box. The administrator password is 123456 by default. Enter the password and click **【OK】** to enter the system setting interface. The user can set the hospital information, device and sample settings, rhythm lead, print, export, display format, display color, filter, FTP (optional), DICOM (optional), Worklist (optional), HL7 (optional), upload (optional), path, other relevant information according to the actual needs.

Name	Function
Hospital information settings	Including hospital title, user-defined term and auto fill .
Device & Sample Settings	<ol style="list-style-type: none"> <li>1) The Device, COM NO., Device NO., Sample Mode, PreSampling Time and Sensitivity of Pace Detection can be set.</li> <li>2) Device: USB ECG device and Bluetooth ECG device. COM No.: select the type of Bluetooth outgoing COM port; only set when the “device” is Bluetooth ECG Device.</li> <li>3) Device NO.: represents the unique number of the equipment, and can be input in English, numbers and special symbols.</li> <li>4) Sampling Time in Auto Mode: it can be entered manually to control the sampling time in automatic record mode; this setting is only for automatic record mode.</li> <li>5) Sample Mode: manual and auto</li> <li>6) PreSampling Time: 3 seconds, 5 seconds and 10 seconds can be selected. It is off by default.</li> <li>7) Sensitivity of Pace Detection: including low, medium, high and off. When the patient has implanted pacemaker, it is necessary to turn on the pacemaker detection sensitivity function, otherwise the pacemaker signal cannot be detected.</li> </ol>

Name	Function
Rhythm leads settings	<p>1) Single rhythm lead: one of the twelve leads is selected as rhythm lead. Single lead rhythm is 1R lead in the format of 3 × 4 + 1R, 6 × 2 + 1R and 1x1.</p> <p>2) Three rhythm leads: three of the twelve leads are selected as rhythm leads. The three rhythm leads are the 3R leads in the format of 3 × 4 + 3R.</p>
Print settings	<p>1) Including horizontal printing report format and vertical printing report format. The horizontal printing format includes horizontal 12x1, horizontal 6x2, horizontal 6x2 + 1R, horizontal 3x4 + 1R, horizontal 3x4 + 3R, and horizontal rhythm lead; the vertical printing format includes vertical 12x1, vertical 6x2, vertical 6x2 + 1R, vertical 3x4 + 1R, vertical 3x4 + 3R, and vertical rhythm lead.</p> <p>2) Print report with ECG grid: enable this item, and the report will be printed with ECG grids.</p> <p>3) Print report with color grid: enable this item, and the report printed will be colorful</p> <p>4) Grid is represented by dots: enable this item and the report will be the grid in the form of dots.</p> <p>5) Enable preview function before print: enable this item to display preview interface when printing report manually.</p> <p>6) Auto-half gain in case of high QRS: if this option is selected, the high gain waveform of chest lead in the printing report will be halved automatically.</p> <p>7) Auto print: enable this item to stop ECG data acquisition and auto print report.</p> <p>8) Thick grid / thin grid: set the color of reporting ECG grid line, which can be modified manually.</p> <p>9) Report title: set the title of the report.</p> <p>10) Enable this item, there will be one more page in the report filled with detailed information of average template.</p>
Export settings	<p>1) Export format: PDF, JPEG, BMP, FDA-XML(optional), EM-XML, SCP (optional), BKG (optional), DICOM.</p> <p>2) Export report display format: horizontal format includes horizontal 12x1, horizontal 6x2, horizontal 6x2 + 1R, horizontal 3x4 + 1R, horizontal 3x4 + 3R, and horizontal rhythm lead; vertical format includes vertical 12x1, vertical 6x2, vertical 6x2 + 1R, vertical 3x4 + 1R, vertical 3x4 + 3R, and vertical rhythm lead.</p>
ECG display format settings	<p>1) Display format: there are six display formats in 12 lead mode: 12x1, 6x2, 6x2 + 1R, 3x4 + 1R, 3x4 + 3R, 1x1.</p> <p>2) Display all paper speed: if this option is checked, all paper speeds will be displayed, and all paper speeds can be switched in the main interface; if this option is not checked, only common paper speeds will be displayed, and only common paper speed can be switched in the main interface.</p> <p>3) Paper speed: set the default paper speed of the ECG waveform in the main interface. All paper speeds include 5 mm / s, 6.25 mm / s, 10 mm / s, 12.5 mm / s, 25 mm / s and 50 mm / S; common paper speeds include 12.5 mm / s, 25 mm / s and 50 mm / s.</p> <p>4) Display all gains: enable to display all gains, display all gains, and all gains can be switched in the main interface ; disable to display all gains, display only common gains, and only common gains can be switched in the main interface.</p> <p>5) Gain: set the default gain of the ECG waveform of the main interface. All the gains include 1.25 mm / mV, 2.5 mm / mV, 5 mm / mV, 10 / 5 mm / mV, 10 / 5 mm / mV, 20 / 10 mm / mV, 20 mm / mV, 20 mm / mV, 40 mm / mV and AGC (automatic gain). Common gains include 5 mm / mV, 10 mm / mV and 20 mm / mV.</p>
ECG display color settings	The background and waveform of the main interface support multiple color selection.
Filter settings	<p>1) AC Filter: provide 50Hz, 60Hz and off 3 options. The purpose of AC filter is to eliminate the interference of AC power supply. The AC filter value should be set according to the natural frequency of the national AC power supply.</p> <p>2) Baseline filter: 7 options are provided: 0.05Hz, 0.16 Hz, 0.25 Hz, 0.32 Hz, 0.5 Hz, 0.67 Hz and off. Baseline filter makes the baseline of ECG waveform regress rapidly, and the larger the filter value, the faster the baseline regression.</p> <p>3) EMG filter: it provides 7 options including 20Hz, 25Hz, 30Hz, 35Hz, 40Hz, 45Hz and off. The smaller the filter value is, the stronger the ability of eliminating EMG interference is.</p> <p>4) Low pass filter: there are four options: 75Hz, 100Hz, 150Hz and off.</p>

Nom	Fonction
ECG Analysis Settings	1) Enable auto-analysis: enable this item, the software can automatically analyze the ECG data collected for more than 10s and provide interpretation. 2) Enable auto-classify: enable this item, the ECG data interpretation of automatic analysis will be classified. 3) Analyze 10s ECG only: enable this item, only the last 10s ECG data collected will be analyzed. 4) QTc duration correction using formula: calculation method of QTc, including two formulas of Bazett and Friderica.
FTP settings (optional)	Includes local FTP server setup and FTP client setup FTP server can receive uploaded file types including EM-XML, FDA-XML, SCP, BKG, DICOM.
DICOM settings (Optional)	According to the corresponding information set by the DICOM server, there are three formats of the uploaded DICOM file which are LEADS_12_ECG, GENERAL_ECG ENCAP-SULETED_PDF_ECG
Worklist settings (Optional)	1) Worklist is used for patient information query: if this function is enabled, the worklist Update button will be displayed in the patient information interface; if this item is not checked, the worklist Update button will not be displayed in the patient information interface. 2) Worklist server setup: set corresponding IP / server port / server aetitle / local AE title information according to worklist server. 3) Worklist update interval in automatic mode: update interval time in automatic update mode, and in manual update mode this button is gray. 4) Update mode: auto, manual
HL7 settings (optional)	1) HL7 is used for patient information query: check this item to open the HL7 query function. Enter the ID in the patient information interface and click <b>【Query】</b> to query the corresponding patient information on the HL7 server according to the ID. If this item is not checked, the HL7 query function will be closed. 2) HL7 server setup: set the corresponding IP / server port according to HL7 server.
Upload settings (optional)	1) Upload mode: auto and manual. 2) Upload method: FTP and DICOM. 3) Upload image file type: PDF, JPEG, BMP, close. If it is selected as close, the image file will not be uploaded. 4) Upload data file format: FDA-XML, EM-XML, SCP, BKG, DICOM, close. Do not upload data files when you choose to close.
Path settings	Set the path to save the software data file.
Miscellaneous settings	1) Enable beep for heart rate: enable this item, the heart beat sound will be generated when ECG waveform is sampled. 2) Lead off alarm: enable this item, when acquire ECG data, there will be a prompt tone when the electrode leads off. 3) Demo mode: check this item, when ECG data is acquired, the main interface will display 84bpm demonstration ECG waveform, and the interface will display the word "demo". In this mode, real-time ECG data will not be acquired. 4) 6 leads mode: check this item, the acquisition interface only displays the limb lead waveform, and the data file and report only save the limb lead data. In this mode, the display format and report format cannot be modified. The display format only supports 6X1 format, and the report only supports horizontal 6X1 format. 5) Delete imported files: check this box to delete the file in the original disk path after import. 6) System password: you can change the administrator password. Click <b>【Change】</b> and enter a new password then click <b>【OK】</b> to complete the password modification.

## EN 5. Instructions for ECG Acquisition and Analysis Software for Android Tablet

### 5.1 ECG Acquisition and Analysis Software for Android Tablet Installation

#### 5.1.1 System Installation and Operation Environment



Android Tablet configuration	Specific Requirement
CPU	Pentium P4
Memory	≥1GB
Storage space	≥4GB
Android version	Bluetooth 4.0 Android ≥4.3
Battery Life	≥6 hours
Screen size	≥7 inch
Resolution	≥1024X600
Printer configuration	Resolution above 600dpi, support network printing function
Others	4.0 Bluetooth

**Note** : Before choosing the android tablet model, please consult the manufacturer or local seller to see if the ECG acquisition and analysis software supports this model.

#### 5.1.2 Software Installation

##### 5.1.2.1 Software Installation

1) Copy ECG acquisition and analysis software for Android Tablet to the root directory of the tablet application device.

2) Open the ECG acquisition and analysis software for Android Tablet, as shown in the following figure, and select «Install»


3) After installation, select **【Open】** to open the license code interface below. After obtaining the license code, enter the license code and click **【OK】** to complete the software license (license prompt will not appear for upgrade installation, please ignore this step).

The screenshot shows a 'Registration' dialog box with two input fields: 'MachineCode' and 'License'. The 'MachineCode' field contains the text '78CA - 0CF0 - 84B4 - C181'. The 'License' field is empty. There are 'OK' and 'Cancel' buttons at the bottom.

## 5.2 Main Interface Function Introduction


### 5.2.1 System Startup







Click the desktop icon to open the ECG acquisition and analysis software for Android Tablet. The initialization main interface is as follows:

### 5.2.2 Introduction of Main Interface and Toolbar

The main interface display includes acquisition and ECG waveform replay. The acquisition interface displays the acquired ECG waveform in real time, as shown below. After the acquisition, enter the replay interface, and move the progress bar to review the acquired ECG waveform in the whole process.

Icon	Button Name	Function Introduction
	New patient	Click <b>【New patient】</b> to open the patient information dialog box, enter the new patient information, and click <b>OK</b> to start sampling ECG data for the new patient.
	Sampling	Click <b>【Sampling】</b> to acquire ECG data directly.
	Stop	Click <b>【Stop】</b> to stop ECG data acquisition.
	Paper speed	Click <b>【Paper speed】</b> to adjust the paper speed of waveform.
	Gain	Click <b>【Gain】</b> to adjust the waveform gain. Note: when the gain is set to 10 / 5 mm / mV, 20 / 10 mm / mV and AGC (automatic gain), the gain of the main interface is displayed as 10 mm / mV. If the gain of the channel is not 10 mm / mV, the multiple relationship is added next to the lead symbol, where (x2), (x4), (x1 / 2), (x1 / 4), (x1 / 8) are respectively represented as the gain of the current channel as 20 mm / mV, 40 mm / mV, 5 mm / mV, 2.5 mm / mV and 1.25 mm/mV.
	Filter	Click <b>【Filter】</b> to adjust the current ECG filter parameter value during acquisition. Filter includes AC filter, EMG filter, baseline filter and low pass filter.
	Setting	Click <b>【Setting】</b> to enter the system setting interface, where you can set related settings.
	Data base	Enter the database, you can query, edit, export, upload (optional) and other operations for the record patient record.
	Exit	Exit the system.

### 5.3 ECG Acquisition and Report Print

The ECG acquisition and analysis software for Android Tablet can acquire ECG data through USB wired acquisition and Bluetooth wireless acquisition. Before sampling ECG data, please make sure that the patient has connected the electrode, the patient cable correctly connected to the ECG acquisition box, and the acquisition box reliably connected to the tablet. See [3.Before Operation] for details.


#### ! Note :

- 1) In case of waveform pause or acquisition interruption during acquisition, please check the reliability of connection between acquisition box and software.
- 2) When the waveform interference is large or the lead falls off during the acquisition, please check the reliability of the connection between the patient cable and the acquisition box and the human body.
- 3) In the acquisition process, if the distance between software and Bluetooth ECG acquisition box exceed the effective transmission distance, the acquisition will be interrupted automatically.

#### 5.3.1 Input Patient Information

Open the ECG acquisition and analysis software for Android Tablet, click **【New Patient】** on the toolbar, enter the corresponding information in the pop-up “Patient Information” dialog box, and click **【OK】** :

The patient information interface includes the following contents:

Name	Function
Patient Information	You can input patient ID, name, age, birthday and gender.
Worklist (Optional)	Enable Worklist in system settings to download patient information from the specified server. Click <b>【Worklist】</b> to enter the Worklist list interface, click <b>【Update】</b> on the list page to enter the update interface, and query the patient information in the server according to the selected examination date. Click patient information in the acquired list to acquire ECG data with this information.
HL7 (Optional)	System setting - hospital information, turn on the HL7 switch, enter the ID in the patient information interface, click <b>【HL7】</b> , and automatically query the patient information of the corresponding ID in the corresponding server. Save the current input patient information and enter the main interface to start acquiring ECG waveform.
OK	Save the current input patient information and enter the main interface to start sampling.
Cancel	Do not save the current input information and return to the main interface.

### 5.3.2 ECG Acquisition

#### 1. Start sampling

Click **【OK】** in the patient interface to start sampling ECG data, or click **【Sampling】** in the toolbar to directly sample ECG data.

Presampling: Set the presampling time of “Settings”- > “Record Mode” to the non close state. ECG data were not recorded in this setup.

After the set pre acquisition time, the system starts to record ECG data, and the recording time of the main interface starts to count.

In the acquisition process, click **【Sampling】** to clear the recording time and restart recording ECG data for the current patient.

It is strongly recommended that when there is strong interference in the current acquisition process, click [Sampling] again to eliminate the interference after the interference disappears, so as to obtain high-quality ECG waveform and reduce misdiagnosis.

#### 2. Stop Sampling

Manual stop: click **【Stop】** to end the current acquisition. When the recording time is less than 10s, the software does not save the data.


Auto stop: if you select «auto» in «Settings» - >»Record Mode» - >»Record Mode», the ECG acquisition and analysis software for Android Tablet will automatically stop and save the file after the record reaches the set sampling time.

## EN 5.3.3 Lead off Alarm

During the acquisition process, if the lead connection is unstable or dropped, the software will prompt the alarm of lead off. Please check the loose electrode or patient cable according to the information prompted, and reconnect the electrode or patient cable.

## 5.3.4 Reanalysis

After the acquisition, if there is a large difference between the measurement parameters provided by the system and the actual waveform, users can click the «Reanalysis» to enter the reanalysis interface, as shown below.

The reanalysis interface shows ECG features, average template of each lead, analysis interpretation, etc., which can be used for re-analysis of current ECG records.

Modify parameters: users can manually input parameters or move the location line of the average template to modify the parameters.

Switch display channel: users can click the lead name to switch to display the average template of the corresponding lead.


Modify conclusion:

- Add: directly input the conclusion in the conclusion input box manually or select the conclusion in the analysis result index, and click [Add] to add the conclusion.
- Clear: clear the conclusion information.
- Save as template: save the current conclusion as user-defined template.
- Edit template:
  - ◆ Import template: import user-defined templates in bulk.
  - ◆ Update: update the selected template content.
  - ◆ Save as template: save the edited template as a new template.
  - ◆ Delete: delete the selected template.
  - ◆ Close: close the edit template box.

## 5.3.5 File Export and Print

### 5.3.5.1 Manual Print and Export

After sampling ECG data, long press ECG waveform replay area to display file export and report print dialog box, as shown below.

Name	Function
Export	Image file format: PDF, JPEG, close. Data file format: EM-XML, FDA-XML (Optional), SCP (Optional), BKG (Optional), DICOM. Image report mode: select the image display format, and the image report format generated after successful export is the selected image report format. Image report preview: automatically open the image report after the export is successful.
Print	Select “report print” to enter the report print interface, select the corresponding report mode, click <b>【OK】</b> , select the corresponding network printer, and click <b>【Print】</b> to print the current report.

There are two ways to implement manual export:

Manual export of replay interface	1) Select “file export” to enter the “file export” interface; 2) Select the exported data file format, image type, and exported image report format. 3) Click <b>【OK】</b> to export.
Database manual export	1) Select one or multiple patient records to export, and click <b>【Export】</b> . 2) In the export file dialog box, select the format of the exported data file. 3) Click <b>【OK】</b> to display the export progress bar. After the file export is completed, the progress bar will be closed and the file export status will be updated.

**! Note:**

- 1) Database export only supports data file export, not image export.
- 2) Database export can realize batch export.


*5.3.5.2 Automatic export and print*

Automatic export: «Setting» → «Record Mode», select «Auto»in export mode , and the report will be exported automatically after the acquisition.

Automatic print: «Settings» → «Report», turn on «auto print», and enter the print preview interface after acquisition.

## EN 5.4 Database

Click **【Database】** in the main interface to enter the database interface. In this interface you can open, query, delete, edit, upload (optional), export patient records. The database interface is shown below.

The database interface functions are introduced in the following table:


Name	Function
Open the patient record	Directly click the patient record to be opened, enter the main interface and display the ECG waveform and patient information of the selected patient record.
Query	Query criteria include ID, name, age, gender, PC synchronization, upload status, export, check date, etc. You can set corresponding query criteria through a single or multiple criteria combination query, and click <b>【Search】</b> to realize the search function.
Patient record management	All the buttons in the interface are gray by default. Select one or more patient records, all the buttons light up. The interface will display the selected number of patient records and the total number of patient records.
Delete	Delete the selected patient record.
Edit	Select the patient record and click <b>【Edit】</b> to edit the patient information.
Upload (optional)	Select one or more patient records and click <b>【Upload】</b> to upload the selected patient records to the corresponding server according to the settings.
Export	Select the patient record to be exported, and click <b>【Export】</b> to open the export file dialog box. Select the format of the exported data file, and click <b>【OK】</b> to display the export progress bar. After the file is exported, the progress bar will be closed and the file export status will be updated.

## 5.5 System Settings

Click **【Setting】** to open the administrator password input dialog box (default password 123456). Enter the password and click **【OK】** to enter the system settings interface. The system setting includes hospital information setting, ECG record, rhythm leads, reports, display format, display colors, filter, and other settings. The user can set the relevant information according to the actual needs.

### 5.5.1 Hospital Information

You can set hospital name, upload settings (optional), Worklist (optional), HL7 (optional) and other relevant information in the hospital information setting interface, as shown in the following figure.


**Note:** HL7 and Worklist are optional functions, and they can not be enabled at the same time. Select one to open and the other to close automatically.

Name	Function
Hospital Information	The name of the hospital is reflected in the report title.
Upload Mode (optional)	Manual: users need to upload data manually. Automatic: sample ECG record and upload it to corresponding server automatically
Upload Method	FTP: select the FTP upload mode, and set the relevant information of FTP server correctly. DICOM: select the DICOM upload mode, and set the relevant information of DICOM server correctly.
Data File Upload Format (optional)	Display only when FTP is selected as upload mode, including EM-XML, FDA-XML, SCP, DICOM and BKG.
Image File Upload Type (optional)	Uploaded image files, including PDF / jpeg.
FTP Server Address/ FTP Server Port/ FTP User Name / FTP Password (optional)	Input corresponding information according to FTP server; when FTP is selected as upload mode, ECG acquisition software will upload the file to FTP server as client.
DICOM Server Address/ DICOM Server Port/ Remote AE Title/ Local AE Title/ DICOM Upload Format (optional)	Input the corresponding information according to the DICOM server; select the DICOM upload mode to upload the patient record data to the corresponding DICOM server.
Worklist Server Address/ Worklist Server Port/ Remote AE Title / Local AE Title (optional)	Set corresponding information according to Worklist server.
Clear Worklist (optional)	Clear the patient information queried from the Worklist server.
Worklist On (optional)	Open this item, and the patient information interface displays the Worklist query function.
HL7 Server Address / HL7 Server Port(optional)	Set corresponding information according to HL7 server.
HL7 On (optional)	Open this item, the patient information interface displays the HL7 query function.

## EN 5.5.2 Record Mode


The machine number, record mode and related exported information can be set in this interface, as shown in the following figure.

Name	Function
Machine Number	Represents the unique number of the software, which can be input in English, numbers and special symbols. The machine number can be reflected in the report name.
Record Mode	Manual: after sampling, the user needs to manually click <b>【Stop】</b> to end the acquisition. Auto: after the start of sampling, when the time reaches the preset sampling time, the software automatically ends the acquisition.
Auto Record Duration	It can be entered manually to control the sampling duration under the automatic recording mode; this setting item is only for the automatic record mode.
Presampling	The pre-sampling time can be 3 seconds, 5 seconds and 10 seconds. It is off by default.
Export Mode	Auto: after ECG data acquisition, the software automatically exports the ECG data file to the corresponding path. Manual: manual data export is required.
Data File Format	Set the default export data file type.
Image File Type	Sets the default export image file type.
Image File Preview	Turn on this item, and the image report will be opened automatically after the image report is exported.

## 5.5.3 Rhythm Lead

Single rhythm lead and three rhythm leads can be set in twelve lead mode in rhythm lead setting interface, as shown below.


Single rhythm lead: one of the twelve leads is selected as rhythm lead. The single rhythm lead is 1R lead in the format of  $3 \times 4 + 1R$  or  $6 \times 2 + 1R$ .

Three rhythm leads: three of the twelve leads are selected as rhythm leads. The three rhythm leads are the 3R leads in the format of 3 × 4 + 3R.

### 5.5.4 Report


In the report interface, you can set the relevant information for printing the report, as shown in the following figure:

Name	Function
Report Format	Set the default report format, including horizontal and vertical report format. The horizontal including horizontal 12x1, horizontal 6x2 + 1R, horizontal 3x4 + 1R, horizontal 3x4 + 1R (three lines), horizontal 3x4 + 3R, horizontal 1R (three lines), horizontal 1R (six lines); vertical format includes vertical 12x1, vertical 6x2 + 1R, vertical 3x4 + 1R (one line), vertical 3x4 + 1R (three lines), vertical 3x4 + 3R, and vertical 1R (Three lines) and vertical 1R (six lines).
Paint Grids	Turn this on and the report will be exported with grid.
Paint Color	Turn this on and the grid will be colorful.
Use Dots instead of Grids	Turn this on and the grid in the report will be replaced by dots.
QRS Auto Halve	Turn this on and the gain of chest lead in the report will be halved automatically.
Auto Print	Turn on this option, you will enter the print report preview interface automatically after acquisition.

### 5.5.5 Display Format


Lead display format and default gain and paper speed values of the main interface can be set in the display format setting interface, as shown in the following figure.

Name	Function
Display Format	Set the default display format of the main interface, including six gears: 12x1, 6x2, 6x2 + 1R, 3x4, 3x4 + 1R, and 3x4 + 3R.
Paper Speed	All paper speeds: 5 mm / s, 6.25 mm / s, 10 mm / s, 12.5 mm / s, 25 mm / s, 50 mm /s. Common paper speeds: 12.5mm/s, 25mm/s, 50mm/s, 12.5 mm / s, 25 mm / s, 50 mm /s.
Speed Select	All and common paper speeds can be chosen.
Gain	All gains: 1.25 mm / mV, 2.5 mm / mV, 5 mm / mV, 10 / 5 mm / mV, 20 / 10 mm / mV, 20 mm / mV, 40 mm / mV, AGC (automatic gain). Common gains: 5mm / mV, 10 mm / mV, 20 mm / mV.
Gain Select	All and common gains can be chosen.


### 5.5.6 Display Color

The background of main interface and the color of wave can be set in the display color setting interface.

### 5.5.7 Filter


Default filter values of the acquisition interface can be set in the filter setting interface including AC filter, EMG filter, Baseline filter, Low pass filter, etc.

Filter	Function Introduction
AC Filter	Three options of 50Hz, 60Hz and off are provided. The purpose of AC filter opening is to eliminate the interference of AC power supply. The selection of AC filter value shall be set according to the natural frequency of national AC power supply.
Baseline Filter	There are 7 options: 0.05 Hz, 0.16 Hz, 0.25 Hz, 0.32 Hz, 0.5 Hz, 0.67 Hz and off. Baseline filter makes the baseline of ECG waveform regress rapidly, and the larger the filtering value, the faster the baseline regression.
EMG Filter	There are 7 options: 20 Hz, 25 Hz, 30 Hz, 35 Hz, 40 Hz, 45 Hz and off. The smaller the filtering value is, the stronger the ability of eliminating EMG interference is.
Low Pass Filter	There are four options: 75 Hz, 100 Hz, 150 Hz and off.


### 5.5.8 Auto-Analysis

Auto- Analysis	Function introduction
Auto Analysis	Enable this item, the measurement parameters and analysis interpretation can be automatically generated after more than or equal to 10 seconds' acquisition.
Diagnostic Classification	Enable this item, ECG at the end of acquisition more than or equal to 10s will be classified according to the analysis interpretation. The software includes normal ECG, borderline ECG, atypical ECG and abnormal ECG.

### 5.5.9 Others

Other settings include demo mode, 6-leads mode, lead off alarm, heart beat sound, pacemaker sensitivity and password modification, as shown in the following table:

Name	Function
Demo Mode	Enable this item, click <b>【Sampling】</b> to display 84bpm demonstration ECG waveform, and the words "demo mode" displayed in the interface. In this mode, real-time ECG data will not be acquired.
6-Leads Mode	Enable this item, the acquisition interface only displays the limb leads waveform, and the data file and report only save the limb lead data.
Lead off Alarm	Enable this item, when the electrode falls off, there will an alarm.
Heart Beat Sound	Enable this item, when sampling, there will a sound.
Pacemaker Sensitivity	Including low, medium and high sensitivity. When the patient equipped with pacemaker, it is necessary to turn on the pacemaker detection, otherwise, the pacemaker signal cannot be detected.
Password Modification	Modify the administrator default password.

## 6. Cleaning, Disinfection and Maintenance

### 6.1 General

Please keep your equipment and its accessories free from dust. To prevent damage to the equipment, be sure to follow the following steps.

- ❗ Before cleaning, disinfecting and maintaining the instrument, be sure to disconnect the power supply (for the equipment with Bluetooth function, please make sure that the battery has been removed), and take off the patient cable.
- ❗ When cleaning, please be careful not to let the liquid penetrate into the instrument. In any case, please do not immerse the instrument and the patient cable into the liquid.
- ❗ Do not clean the unit and accessories with abrasive fabric and avoid scratching the electrodes.
- ❗ After cleaning, avoid any detergent residue on the surface of the instrument and patient cable.
- ❗ This instrument and all its accessories must be regularly maintained and overhauled (at least once a year).
- ❗ ECG machines are classified into Measuring Instrument, so users should send it to Official Measurement Administrative to test and certify according to national metrological calibration regulation of electrocardiograph and electroencephalogram machine every year.
- ❗ The instrument signal input/output connector (when needed) must be connected with Class I equipment which is IEC60601-1 compliant, and the total leakage current should be tested to be available by users themselves.
- ❗ This instrument needs to be maintained and repaired by trained and qualified professional personnel; when the instrument functions abnormally, it should be clearly identified to avoid the instrument running with faults.
- ❗ The instrument shall not be remodeled or modified in any way.

### 6.2 Cleaning

#### 6.2.1 Cleaning of Acquisition Box

Wipe the surface of the acquisition box with a cloth dipped in mild neutral detergent (or 75% alcohol). Then wipe it with a dry cloth. Avoid water entering the acquisition box, which may cause malfunction.

#### 6.2.2 Cleaning of Patient Cable and Electrodes

Clean the patient cable and electrodes with a soft cloth dipped in water, neutral detergent or alcohol (75%), and dry them with a dry cloth to ensure that the patient cable and electrodes are completely dry after cleaning.

- ❗ **Note:** if there is gel on the surface of patient cable or electrodes or electrode pins or metal parts are wet, inaccurate recording of ECG waveform may occur.

### 6.3 Disinfection

It is recommended that disinfection should be only performed under the regulation of medical institutions or under other proper situations; Disinfecting equipment must be disinfected before cleaning the ECG machine.

- ❗ Do not use high-temperature, high-pressure steam, ionizing radiation methods for disinfection.
- ❗ Do not use chlorine-containing disinfectant for disinfection, such as chlorinated lime and sodium hypochlorite etc.

After the power is turned off, use a soft and clean lint free cloth to absorb a proper amount of disinfectant (75% alcohol), wipe the external surface of accessories such as ECG acquisition box and patient cable, chest electrode suction ball and limb electrode clamp, and use a dry cloth to wipe off the residual disinfectant. After disinfection, please place the ECG acquisition box and its accessories in a cool and ventilated place to dry thoroughly.

### 6.4 Care and Maintenance

- ❗ This instrument needs to be maintained and repaired by trained and qualified professional personnel; when the instrument functions abnormally, it should be clearly identified to avoid the instrument running with faults.

#### 6.4.1 General Maintenance

In addition to the care and maintenance methods recommended in this manual, other measurement requirements and maintenance methods specified in local regulations shall be followed at the same time.

Periodic safety test shall be carried out for the instrument. The test cycle shall be at least once a year. The test contents mainly include:

- a) Check whether there is mechanical and functional damage in the main unit and accessories;
- b) Check whether the safety signs are damaged;
- c) Verify the function of the instrument according to the Instructions for Use;
- d) The following safety tests are carried out in accordance with EN 60601-1:2006:
  - Earth leakage current, limit: NC 50 $\mu$ A, SFC 1000 $\mu$ A
  - Enclosure leakage current, limit: NC 100 $\mu$ A , SFC 500 $\mu$ A
  - Patient leakage current, limit: 10  $\mu$ A (CF type equipment)
  - Patient auxiliary current, limit: NC A.C. 10 $\mu$ A, D.C. 10 $\mu$ A; SFC A.C. 50 $\mu$ A, D.C. 50 $\mu$ A
  - Patient leakage current under single fault condition with mains voltage on the applied part, Limit: 50  $\mu$ A (CF Type)

Inspection and record shall be conducted by professional personnel with safety test knowledge and qualified training, and practical experience to perform these tests. If the instrument has any problems in the above tests, it must be repaired.

#### 6.4.2 Care and Maintenance of ECG Acquisition Box

- Take out the battery to make the acquisition box in power off state (for PE-1204 and PE-1204A);
- Clean the instrument, wipe the accessories and cover the dust cover;
- Prevent liquid from penetrating into the interior of the instrument, affecting the performance and safety of the instrument;
- Ask the medical instrument maintenance department to verify the instrument regularly.

### 6.4.3 Care and Maintenance of Patient Cable

- Check the patient cable is well-collected according to the following table. The resistor between the electrode plug and plug pin of the patient cable plug should be less than  $10\ \Omega$ .
- Note: the resistance or of patient cable with the function of defibrillation is around  $10K\ \Omega$ .

Wire PlugSymbol	R	L	F	RF	C1	C2	C3	C4	C5	C6
Plug Pin Location	9	10	11	14	12	1	2	3	4	5

- Check the integrity of patient cable;
- Knotting or crooking in closed angle will shorten the patient cable's life, so please align the patient cable and then collect the electrodes;
- Do not drag or twist the patient cable with excessive stress while using. Hold the connector plugs instead of the cable when connect or disconnect the patient cable;
- Storing the lead wires in bigger wheel or swinging to avoid twisting or folding;
- Once damaged or aged patient cable has been found, replace it with new one immediately.

### 6.4.4 Care and Maintenance of Electrodes

- Electrodes must be cleaned after use and be sure there is no remaining gel on them.
- Keep the suction ball of chest electrodes away from sunshine and excessive temperature.
- After long time use, the surface of electrodes will be oxidized because of erosion and other causes. In this case, electrodes should be replaced to achieve high-quality ECG.

## 7. Common Troubleshooting

### 7.1 Common Troubleshooting of ECG Acquisition and Analysis Software for Windows PC

Number	Problem	Solution
1	Double click the icon of ECGApp software on the desktop of the computer but cannot access the ECG acquisition and analysis software for Windows PC	<ol style="list-style-type: none"> <li>1. The software in the installation directory may be deleted by mistake. Re-install ECG acquisition and analysis software for Windows PC.</li> <li>2. The computer has no Dongle. Insert the Dongle.</li> <li>3. The USB interface inserted into the Dongle is damaged. Change another USB port and insert the Dongle.</li> <li>4. If the USB port has been changed, the software still cannot be opened ,then change another Dongle.</li> </ol>
2	Communication failure between computer and acquisition box	<ol style="list-style-type: none"> <li>1. Check whether the acquisition box driver is installed successfully, and open the device manager to check whether the acquisition box port displays a yellow exclamation mark. If the yellow exclamation mark is displayed in front of the port of the acquisition box, please install the acquisition box driver manually. For the installation method, please refer to number 16 of common troubleshooting of ECG acquisition and analysis software for Windows PC</li> <li>2. If the USB port of the computer is damaged, change another USB port to insert the data cable.</li> <li>3. If the data cable is damaged, change it.</li> <li>4. If the acquisition box is damaged, change it.</li> </ol>

Number	Problem	Solution
3	Communication failure between computer and Bluetooth ECG acquisition device	<ol style="list-style-type: none"> <li>1. Check whether the Bluetooth ECG acquisition box is opened normally.</li> <li>2. Check whether the Bluetooth ECG acquisition box and the computer are matched successfully.</li> <li>3. Check whether the COM number of the device is set correctly, and whether the COM number of the device is consistent with the Bluetooth outgoing COM port.</li> <li>4. Check whether the battery power of Bluetooth ECG acquisition box is sufficient. If not, replace it with a new battery.</li> <li>5. Shorten the transmission distance between Bluetooth ECG acquisition box and computer.</li> <li>6. If the Bluetooth ECG acquisition box is damaged, change it.</li> </ol>
4	What should you do when Bluetooth transmission acquisition is interrupted suddenly or the signal is unstable?	<ol style="list-style-type: none"> <li>1. Check whether the battery power is sufficient; change the battery when the battery power is insufficient.</li> <li>2. Ensure that the surrounding environment is free of large electromagnetic interference and other interference of Bluetooth devices.</li> <li>3. Shorten the transmission distance between Bluetooth ECG acquisition box and computer.</li> </ol>
5	After pressing 【Print】 , the printer did not respond after a few minutes.	<ol style="list-style-type: none"> <li>1. Whether the default printer in the computer is set.</li> <li>2. If the default printer is set, check whether the printer power is turned on and whether the USB data cable of the printer is connected normally.</li> <li>3. If the power supply and data cable are connected normally. Check the printer status on the right side of the start menu to see if there are too many print jobs waiting in the printer status bar.</li> <li>4. If the printing job is piled up and cannot be printed, it is recommended to restart the computer and print again.</li> <li>5. If still fail to print normally, the printing button is pressed after the equipment is in normal use for a period of time and the software reports an error and has no response, it may be that some security software installed on the customer's computer mistakenly deleted the files in the software. Re-install the software again.</li> </ol>
6	File import failure	Check whether the imported data file is the ECG data file format (EM-XML, FDA-XML, SCP, DICOM, BKG) generated by the software provided by ECGMAC.
7	The exported image report is not clear and blank, or there is a problem in the display	Check whether the default printer in the computer is set or not. Set the normal printer or Microsoft XPS printer as the default printer in the computer.
8	In the process of ECG acquisition, interference often appears in some fixed leads in 12 lead ECG..	<ol style="list-style-type: none"> <li>1. Check the patient cable, unplug it from the acquisition box and reconnect it again</li> <li>2. Check whether the electrode is in poor contact. The human body part contacted by the electrode needs to be wiped with alcohol and then reconnected.</li> <li>3. If the problem still can't be solved, check whether the suction ball or clamp corresponding to the disturbed lead is rusted at the metal connection. After treatment, reconnect them again</li> <li>4. If the problem cannot be solved, change the patient cable.</li> </ol>

Number	Problem	Solution
9	In the process of ECG acquisition, ECG baseline wander is serious, and there is burr like interference on the baseline.	<ol style="list-style-type: none"> <li>1. First, modify the filter parameter value in the software. Common filter parameter value of baseline filter is 0.67 Hz, EMG filter 40 Hz, AC filter 50 Hz.</li> <li>2. Check the patient cable. If the appearance of the patient cable is normal, change it and then conduct sampling test. (Note: the patient cable of the acquisition box and the patient cable of other electrocardiographs are the same, both of which are DB25 pin interfaces. Find a patient cable that can be used normally to replace the old one, and check whether it is the problem of ECG lead).</li> </ol>
10	The operating system is too slow or has other problems. Therefore, User wants to reinstall the operating system. However, the patient record information acquired before is also needed to be preserved for later check. How to solve it?	<p>Before re-installing the operating system, if you want to save the previous patient records, copy the data file to a safe storage device. After the system re-installation, enter the patient record, click <b>【Import】</b> and select the backup data file.</p>
11	Aucune forme d'onde après la connexion du câble patient et du corps humain	<ol style="list-style-type: none"> <li>1. Check whether there is limb lead or chest lead is not connected (you can check the software's lead off information).</li> <li>2. If the software prompts that all leads fall off, please check whether the RF (n) of the right lower limb lead is well connected.</li> <li>3. The electrode is not coated with conductive medium (water, alcohol).</li> <li>4. The patient cable is damaged. Please change another one.</li> <li>5. After the above faults are eliminated, there are still the same problems. The equipment may be damaged. Please return to the factory for maintenance.</li> </ol>
12	The displayed heart rate is twice the actual heart rate of the waveform	<p>Please check whether the connection between the software and the device is matched. The software device is USB ECG device. While the default COM number is consistent with the outgoing port of the Bluetooth ECG acquisition box that has been successfully paired with the computer, USB acquisition can find the Bluetooth acquisition box with the same COM number. However, the waveform will be distorted. The shortening of the R-R interval will cause the displayed heart rate value to be twice the actual heart rate of the waveform.</p>
13	Patient record upload failed, please check FTP settings	<ol style="list-style-type: none"> <li>1. Check if the computer is connected to the Internet.</li> <li>2. Setting - whether FTP upload mode is selected for upload setting.</li> <li>3. Check whether the setting FTP setting account / password / server IP address / port is consistent with the FTP server.</li> <li>4. If SSL encryption upload is selected, you should check whether the selected SSL certificate is valid.</li> <li>5. Check whether the folder receiving data of FTP server exists and whether the read-write permission of the folder is enabled.</li> <li>6. Check whether FTP server is blocked by firewall.</li> </ol>

Numéro	Problème	Solution
14	PC-ECG failed to receive file as local FTP server	<ol style="list-style-type: none"> <li>1. Check whether the FTP server parameters set by the client are correct.</li> <li>2. Check whether the folder receiving data exists.</li> <li>3. Check whether the ECG file service is running normally.</li> <li>4. Run ECGApp as administrator and keep the software on during the upload process.</li> <li>5. Close the firewall or add the ECG file service as the application allowed by the firewall.</li> </ol>
15	How to set ECG file server as the allowed application of firewall?	Right click the computer icon on the desktop → properties → control panel → system and security → windows firewall → allow applications or functions to pass through windows firewall → allow other applications → browse → find the folder where ECG file server is located and select this program → open and add.
16	How to install the acquisition box driver manually?	<ol style="list-style-type: none"> <li>1. Please turn off anti-virus software or computer Butler and other similar software before installing the driver.</li> <li>2. Connect the acquisition box to the USB port of the computer through the USB cable. Windows will prompt to find the new hardware and enter the new hardware installation wizard to search for the driver of the ECG acquisition box.</li> </ol>

## 7.2 Common Troubleshooting of ECG Acquisition and Analysis Software for Android Tablet

Number	Problem	Solution
1	Double click the icon of App on the desktop of the computer but cannot access the software.	<ol style="list-style-type: none"> <li>1. Check the Android system version. The software may not be compatible with the latest system.</li> <li>2. Please check whether all authorities in the application authorities of ECG Acquisition and Analysis Software for Android Tablet are enabled.</li> </ol>
2	Communication failure between tablet and USB ECG acquisition box	<ol style="list-style-type: none"> <li>1. Check whether the connection between USB cable and tablet is loose, unplug and reconnect again.</li> <li>2. If the data line or OTG line is loose, reconnect it.</li> <li>3. If the data cable or OTG line is damaged, replace the data cable or OTG and reconnect again.</li> <li>4. If the acquisition box is damaged, change it.</li> <li>5. If the USB port of the tablet is damaged, change the tablet.</li> </ol>
3	Computer (or tablet) failed to match with Bluetooth ECG acquisition device	<ol style="list-style-type: none"> <li>1. Check whether the Bluetooth ECG acquisition box is opened normally.</li> <li>2. Check whether the battery power of Bluetooth ECG collection box is sufficient. If the battery power is insufficient, replace it with a new battery.</li> <li>3. Shorten the transmission distance between Bluetooth ECG acquisition box and computer (or tablet).</li> <li>4. Check whether the Bluetooth ECG acquisition box is damaged. If it is damaged, change it.</li> <li>5. Change the computer (or tablet).</li> </ol>
4	What should you do when Bluetooth transmission acquisition is interrupted suddenly or the signal is unstable?	<ol style="list-style-type: none"> <li>1. Ensure that the surrounding environment is free of large electromagnetic interference and other interference of Bluetooth devices.</li> <li>2. Replace the battery of the acquisition box.</li> <li>3. Shorten the transmission distance between Bluetooth ECG acquisition box and computer (or tablet).</li> </ol>

Numéro	Problème	Solution
5	The exported image report cannot be opened	If PDF reader is installed on the tablet, re-install PDF reader.
6	Patient record upload failed, please check FTP settings	<ol style="list-style-type: none"> <li>1. Check if the tablet is connected to WiFi.</li> <li>2. Check whether the account / password / server IP address / port set in the selected upload method is consistent with the server.</li> </ol>
7	In the process of ECG acquisition, interference often appears in some fixed lead of 12 leads.	<ol style="list-style-type: none"> <li>1. Check the patient cable, unplug the patient cable from the acquisition box and reconnect it again.</li> <li>2. Check whether the electrode is in poor contact. The human body part contacted by the electrode needs to be wiped with alcohol and then reconnected.</li> <li>3. If the problem still can't be solved, check whether the suction ball or clamp corresponding to the disturbed fixed lead is rusted at the metal connection. After treatment, reconnect again.</li> <li>4. If the problem cannot be solved, replace the patient cable.</li> </ol>
8	In the process of ECG acquisition, ECG baseline wander is serious, and there is burr like interference on the baseline.	<ol style="list-style-type: none"> <li>1. Modify the filter parameter value in the software, common filter parameter value of baseline filter is 0.67hz, EMG filter 40Hz, AC filter 50Hz.</li> <li>2. Check the patient cable. If the appearance of the lead wire is normal, replace the lead wire and then conduct sampling test. (Note: the patient cable of ECG acquisition box and other patient cable of ECG machine are the same, both of which are DB25 pin interfaces. Replace it with another one to check whether it's the problem of patient cable).</li> </ol>
9	No waveform after the connection between patient cable and human body	<ol style="list-style-type: none"> <li>1. Check whether the lead ECG acquisition box is connected with the patient cable.</li> <li>2. Check whether there is any lead not connected (you can check the software's lead off alarm).</li> <li>3. If the software prompts that all leads fall off, please check whether the RF (n) of the right lower limb lead is well connected.</li> <li>4. The electrode is not coated with conductive medium (water, alcohol).</li> <li>5. The lead wire is damaged and should be replaced.</li> <li>6. If there are still the same problems after the above troubleshooting, the equipment may be damaged. Please return to the factory for maintenance.</li> </ol>
10	Network printer can not print report	<ol style="list-style-type: none"> <li>1. Check whether the printing service software of the network printer is correctly installed on the tablet, and refer to the instructions of the printer.</li> <li>2. Whether the power of the network printer and WiFi are turned on.</li> <li>3. Whether WiFi connection is turned on for the tablet.</li> <li>4. Whether the tablet and network printer are in the same LAN.</li> <li>5. Whether the printing service of network printer is enabled on the tablet.</li> </ol>

## 8. Warranty and After-sale Service

### 8.1 Warranty

If a product covered by this warranty is determined to be defective because of defective materials, components, or workmanship, and the warranty claim is made within the warranty period, ECGMAC will, at its discretion, repair or replace the defective part(s) free of charge.

- Main Unit:

ECGMAC warrants that ECGMAC's products meet the labeled specifications of the products and will be free from defects in materials and workmanship occur within warranty period. Within one year from the date of delivery, the manufacturer shall provide a free maintenance after receiving Failure Proof Report.

- Accessories:

Under normal use and maintenance condition, the warranty will be free if a failure report is received within half a year from the date of delivery.

- Software:

The warranty will be free if a failure report is received within a year from the date of delivery. If you need to upgrade the software, please consult the manufacturer directly.

Note: ECGMAC's obligations under this guarantee do not include freight and other costs. All costs for products beyond the warranty period and for product maintenance shall be borne by the user. All maintenance shall be carried out by engineering and technical personnel approved by ECGMAC.

### 8.2 Responsibility of Exemption

ECGMAC is not responsible for the damage and delay caused by the following circumstances:

- The user disassembles, stretches and readjusts the instrument.
- Damage caused by alteration or repair by anyone not authorized by ECGMAC.
- Replacement or removal of serial number label and manufacture label.
- Damage caused by mishandling.
- Damage caused by abnormal use beyond the specified use conditions.

### 8.3 After-sale Service

If the user has any questions, please contact your local distributor or manufacturer.

After-sales service agent: Shenzhen ECGMAC Medical Electronics Co., Ltd.

Residence: 2nd Floor of Block 2, Haoye Industrial Park, Tiegang Road, Xixiang Street, Baoan District, 518102 Shenzhen, China.

Postal Code: 518102

Telephone: +86-755-27697821 0755-27697823

Fax: 0755-27697823

E-mail: info@ecgmac.com

### 8.4 Service Life and Production Date

Production date: see product label for details.

Service life of the product: 5 years (the service life is only limited to the main unit excluding the replaceable accessories; if the product needs to continue to be used after more than 5 years, it needs to be maintained every six months).

## Appendix A: Package and Accessories

### A.1 Accessories

When the product is dispatched from the factory, completed package should contain the following accessories:

PE-1201/PE-1201A/PE-1204/PE-1204A Acquisition Box	1
Software Installation CD	1
USB DONGLE (License Code)	1
Patient Cable	1
Limb electrodes (clamp)	4/set
Chest electrodes (suction ball)	6/set
USB data cable	1
User manual	1

### A.2 Caution

- 1) Open the box from the top;
- 2) After opening the box, check the accessories and user manual, and then check the instrument.;
- 3) If the package does not match the packaging list or the instrument does not work properly, contact the sales department or customer service department;
- 4) Please use the accessories supplied by ECGMAC. Accessories from other suppliers may damage the instrument and affect its performance and safety. Before using the accessories from other suppliers, please consult our customer service first;
- 5) To enable us to serve you in time, please fill out the warranty card (copy) and mail it to us;
- 6) Please keep the instrument package for periodic testing or equipment maintenance.

## Appendix B: Technical Specifications

### B.1 Performance Specifications

Power supply	PE-1201/PE-1201A :d.c 5.0 (power supply by computer USB port )
	PE-1204/PE-1204A :d.c5.0V(power supply by computer USB port )/d.c3.0V ( AAA Alkaline batteries)
Lead	Standard 12 Leads
Input method	Floating Ground input with defibrillator protection
Calibration voltage	1mV±2%
AD converter	12bit( PE-1201/PE-1204)
	24bit( PE-1201A/PE-1204A)
Input circuit current	≤0,01 μ A
	≤0,05
Sensitivity	10 mm/mV ± 2%
Noise level	≤15 μ V (PE-1201/PE-1204)
	≤10 μ V (PE-1201A/PE-1204A)
Polarization voltage	± 400 mV (PE-1201/PE-1204)
	± 1180 mV (PE-1201A/PE-1204A)
Input impedance	≥100MΩ (10Hz) (PE-1201A/PE-1204A)
	≥50MΩ (10Hz) (PE-1201/PE-1204)
Input voltage range	PE-1201/PE-1204: each lead ≥(-7.5-7.5)mV. No distortion of waveform
	PE-1201A/PE-1204A: each lead ≥(-22.5-22.5)mV. No distortion of waveform
CMRR	PE-1201A/PE-1204A : ≥110 dB (AC filter off) ; ≥120 dB (AC filter on)
	PE-1201/PE-1204: ≥100 dB (AC filter off) ; ≥120 dB (AC filter on)
Frequency response	0,05 Hz~150 Hz (PE-1201/PE-1204)
	0,01 Hz~350 Hz (PE-1201A/PE-1204A)
Time constant	≥3,2 s (PE-1201/PE-1204)
	≥5 s (PE-1201A/PE-1204A)
Paper Speed	5 mm/s、 6.25 mm/s、 10 mm/s、 12.5 mm/s、 25 mm/s、 50 mm/s ± 2 %
HR range	HR range: 30bpm~300bpm, with error of ± 1 or ± 1% of the displayed value, whichever is greater
Gain	1,25 mm/mV、 2,5 mm/mV、 5 mm/mV、 10 mm/mV、 20 mm/mV、 40 mm/mV、 10/5 mm/mV、 20/10 mm/mV、 AGC ± 2 %
Sampling rate	2000/Sec/channel (PE-1201/PE-1204)
	32000/Sec/channel (PE-1201A/PE-1204A)
Amplitude quantisation	4,563 uV/LSB (PE-1201/PE-1204)
	2,289 uV/LSB (PE-1201A/PE-1204)
Bluetooth transmission distance	The Bluetooth limited transmission range of Battery powered model ≤ 5m
Operation life of battery	For battery powered models, the continuous working time under the condition of internal battery power supply shall not be less than 240Min.
Pacemaker pulse detection	±2mV~±700mV , Time Duration: 0.5 ms~2.0 ms (PE-1201/PE-1204)
	±2mV~±700mV , Time Duration: 0.1 ms~2.0 ms (PE-1201A/PE-1204A)

## EN B.2 Function Features

ECG Input channel	Simultaneous acquisition of standard 12 leads ECG
Sample mode	6 and 12 lead acquisition mode
Sample method	Support USB and Bluetooth acquisition mode (PE-1201/PE-1201A only support USB acquisition)
AVL, aVR leads correction	Support the correction of left and right upper limb leads to reduce the secondary acquisition caused by misoperation (only ECG acquisition system for Windows system)
Replay	Support data replay
Pacemaker detection	Support pacemaker detection and the sensitivity of the pacemaker detection can be adjusted.
High precision measuring tools	Support ECG scale measurement, parallel ruler, magnifier (only ECG acquisition system for Windows system)
Working mode	Two modes: manual and auto
Rhythm	At least 300 second rhythm acquisition function (only ECG acquisition system for Windows system)
Filter	Four kinds of digital filters which can filter out all kinds of interferences and ensure ECG waveform accuracy a) AC filter: 50Hz, 60Hz, off; b) Baseline filter: 0.05 Hz, 0.16 Hz, 0.25 Hz, 0.32 Hz, 0.5 Hz, 0.67 Hz, off; c) EMG filter: 20Hz, 25Hz, 30Hz, 35Hz, 40Hz, 45Hz, off d) Low pass filter: 75 Hz, 100 Hz, 150 Hz, off.
ECG lead display	Seven display modes: 3×4, 3×4+1R, 3×4+3R, 6×2, 6×2+1R, 12×1, 1R
Sample mode	Presampling and real-time sampling
Data mode	ECG mode and Demo mode
Input method	Floating Ground input with defibrillator protection
Lead system	Standard 12 lead and 6 lead
Defibrillation protection	With anti defibrillation and shock protection function
Patient	Manually input new patient information or query patient information through Worklist ( optional) / HL7 (optional) protocol and other networks
Patient record management (ECG acquisition and analysis software for Windows PC)	A. It has the functions of editing, querying, listing, opening, deleting, comparing, uploading (optional) and sampling a) Edit: edit new patient information; b) Query: query according to query criteria such as ID, name, age, gender, birthday, Department, examining doctor, examining doctor, inpatient number, outpatient number, examination date, number of leads, user-defined information items, etc; c) List: display the inspection time, ID, name, age, gender, uploaded, number of leads, file path, view status and other information of the patient record in the form of list; d) Open: open the selected patient record, and display the ECG data and patient information of the patient record; e) Delete: delete the authority of selecting patient record; f) Comparison: display two patient record comparison functions in the same screen g) Upload: upload the selected patient record (optional) h) Acquisition: use the existing patient information to enter the collection interface
Data transmission and networking functions	1. It shall be able to access the third party HIS or EMR system through FTP (optional), DICOM (optional) and other general standard protocols. 2. It shall be able to have the function of Worklist (optional) and HL7 (optional) to actively obtain the patient application examination information in HIS or EMR system. 3. ECG report and ECG data remote transmission function.

<p>Patient record management (ECG acquisition and analysis software for Windows PC)</p>	<p>B.ECG data import / export function  a) Support the import of ECG data from disk into windows ECG acquisition system  b) Support Android ECG acquisition system to import data into windows ECG acquisition system through USB data cable;  c) Support ECG data to be imported into windows ECG acquisition system through network;  d) Import data format: EM-XML, FDA-XML (optional), SCP (optional), BKG (optional) and DICOM format.  e) Support the export of ECG data from windows ECG acquisition system to disk;  f) Export method: automatic and manual  g) Export image file types: PDF, JPEG, BMP;  h) Export data format: EM-XML, FDA-XML (optional), SCP (optional), BKG (optional) and DICOM format.  C. ECG data upload function (optional):  a) Upload mode: FTP and DICOM  b) Upload method: manual and automatic  c) Upload image file types: PDF, BMP and JPEG  d) Upload data file format: EM-XML, FDA-XML (optional), SCP (optional), BKG (optional) and DICOM (leads_12_ECG, GENERAL_ECG, ENCAPSULETED_PDF_ECG)  e) Upload image display format: 12 lead modes have 12 report formats: horizontal 12x1, horizontal 6x2, horizontal 6x2 + 1R, horizontal 3x4 + 1R, horizontal 1R, vertical 12x1, vertical 6x2, vertical 6x2 + 1R, vertical 3x4 + 1R, vertical 3x4 + 3R, and vertical 1R</p>
<p>Database (ECG acquisition and analysis software for Android Tablet)</p>	<p>a. Data management functions: query, edit, open, delete, upload (optional), export;  b. ECG data export function:  a) Export image file type: PDF, JPEG;  b) Export data format: EM-XML, FDA-XML(optional), SCP (optional), BKG (optional) and DICOM format.  c. ECG data upload function (optional):  a) Upload mode: FTP and DICOM;  b) Upload method: automatic and manual  c) Upload image file type: PDF and JPEG;  d) Upload data file format: EM-XML, FDA-XML (optional), SCP (optional), BKG (optional) and DICOM (DICOM_12L, DICOM_STANDARD, DICOM_PDF).  d. Display format of uploaded image: 12 lead mode has 14 report format including horizontal 12x1, horizontal 6x2 + 1R, horizontal 3x4 + 1R (one line), horizontal 3x4 + 1R (three lines), horizontal 3x4 + 3R, horizontal 1R (three lines), horizontal 1R (six lines), vertical 12x1, vertical 6x2 + 1R, vertical 3x4 + 1R (one line), vertical 3x4 + 1R (three lines), vertical 3x4 + 3R, vertical 1R (three lines) and vertical 1R (six lines);  e. Directly import of data into windows ECG acquisition system through Android ECG acquisition system is supported.</p>
<p>Report preview and report function (ECG acquisition and analysis software for Windows PC)</p>	<p>a. Report format: the 12 lead model has 12 report formats: horizontal 12x1, horizontal 6x2, horizontal 6x2 + 1R, horizontal 3x4 + 1R, horizontal 1R, vertical 12x1, vertical 6x2, vertical 6x2 + 1R, vertical 3x4 + 1R, vertical 3x4 + 3R, and vertical 1R;  b. Auto-half gain in case of high QRS  c. Automatic print is supported  d. AGC is supported;  e. Print preview is supported.  f.User-defined report title is supported</p>
<p>System settings</p>	<p>a. User-defined item of single rhythm lead and three rhythm lead is supported  b. ECG waveform can have multiple colors  c. Setting hospital name is supported  d. Support the setting of heart beat, lead off alarm, etc</p>

Report preview and report function (ECG acquisition and analysis software for Android Tablet)	<p>a. Report format: the 12 lead model has 14 report formats: horizontal 12x1, horizontal 6x2 + 1R, horizontal 3x4 + 1R (one row), horizontal 3x4 + 1R (three rows), horizontal 3x4 + 3R, horizontal 1R (three rows), horizontal 1R (six rows), vertical 12x1, vertical 6x2 + 1R, vertical 3x4 + 1R (one row), vertical 3x4 + 1R (three rows), vertical 3x4 + 3R, vertical 1R (six rows);</p> <p>ECG report output format: ECG grid, colorful report, ECG dot and high QRS waveform gain are automatically halved;</p> <p>b. rhythm lead report mode is supported</p> <p>c. AGC is supported</p> <p>d. automatic preview report after data acquisition is supported</p> <p>e. network printer is supported</p> <p>f. automatic print preview after data acquisition is supported</p>
--	--

### B.3 Safety Standard

EN 60601-1:2006	Medical electrical equipment - Part1: General requirements for basic safety and essential performance
IEC 60601-2-25:2011	Medical electrical equipment - Part 2-25: Particular requirements for the basic safety and essential performance of electrocardiographs
IEC 60601-2-51:2003	Medical electrical equipment - Part 2-51: Particular requirements for safety, including essential performance, of recording and analysing single channel and multichannel electrocardiographs
IEC 60601-1-2	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests
Shell case protection class	IPX0, without waterproof function
Anti-electric-shock type	Special power supply equipment, internal power supply equipment; (when it is connected to the power supply network, the equipment providing power to it shall meet the requirements of class I or class II equipment)
Anti-electric-shock degree:	Type CF with defibrillation-proof
Degree of safety of application in the presence of flammable gas:	Equipment not suitable for use in the presence of flammable gas
Working mode:	Continuous operation

### B.4 Dimension and Weight

Main Unit Dimension	68 mm×126 mm×24 mm
Package Size	70 mm×330 mm×85 mm
Net weight/ Gross weight	0.85 kg /1.7 kg

### B.5 Power Supply Requirements

- 1) USB power supply: Calibration Voltage=5V
- 2) Input Power: 2x1.5V Excell Alkaline AAA Calibration Voltage= 3V

### B.6 Environment Conditions

Transportation	Temperature	Relative humidity	Atmospheric pressure
Storage	-20 °C to +55 °C	≤ 93%	50kPa~106kPa
Operation	-20 °C to +55 °C	≤ 93%	50kPa~106kPa
Fonctionnement	+5 °C to +40 °C	≤80%, sans condensation	86kPa~106kPa

## Appendix C: Key Components

Number	Name	Type	Specification	Note
1	Isolation transformer	EE16	DIP10	
2	Patient cable	ECG-FD08X4	Φ 4/TPU/4-lead 6 channels	
3	Optocoupler	H11L1SR2M or H11L1SM (the suffix represents the different number of packages)	Transmission data rate 1M	

## Appendix D: List of Functional Accessories

Number	Name	Specification	Note
1	Patient cable	ECG-FD08X4	Standard configuration
2	Limb electrode	ECG-FJX41	Optional
3	Chest electrode	ECG-FQX41	Optional
4	Limb electrode	ZJ-01	Optional
5	Chest electrode	XQ-01	Optional
6	Limb electrode	ZJ-02	Optional
7	Chest electrode	XQ-02	Optional
8	OTG data cable	Y-C438	Optional
9	USB cable	1.5m	Standard configuration
10	USB Dongle		Standard configuration

## Appendix E: Network Security Specifications

### E.1 User access control

#### (1) User identification method

User Type	User Identification Method
Administrators	Administrator Password
General Operator	No Administrator Password

#### (2) User type and authority

User Type	User Authority
Administrators	With administrator password, you can enter the setting interface, patient record interface or database interface for setting operation, and have all operation permissions
General Operator	Without administrator password, you only have partial operation permission and you are unable to use setting and operation permission of patient record or database

## E.2 Data exchange mode

Data Exchange type	Data Communication Type	Transfer Protocol	Storage Format
Data exchange between ECG acquisition box and ECG acquisition and analysis software	USB	USB2.0	/
	Bluetooth	Bluetooth4.0	/
Data exchange between ECG acquisition and analysis software for Android Tablet and ECG acquisition and analysis software for Windows PC	USB	USB2.0	EM-XML, FDA-XML, SCP, BKG and DICOM
	Ethernet	FTP	EM-XML, FDA-XML, SCP, BKG and DICOM
Data exchange between ECG acquisition and analysis software and external data	USB/SD	USB2.0	EM-XML, FDA-XML, SCP, BKG, DICOM, BMP, JPG, PDF
	Ethernet	FTP/worklist/dicom/HL7 (Optional)	EM-XML, FDA-XML, SCP, BKG, DICOM, BMP, JPG, PDF (for FTP or DICOM)

## E.3 Operation Environment

Item	ECG Acquisition System for Windows System	ECG Acquisition System for Android System
Hardware Configuration	CPU: $\geq 1.2$ G Memory: $\geq 512$ MB Main board: Inter main board recommended Hard disk: $\geq 40$ GB Display card: $> 1$ M Display: $\geq 43$ cm (17Inch), $> 1024 \times 768$ , $\geq 16$ CD-ROM: $\geq 24$ times Operation system: Microsoft Windows 2000/xp/7/8/10 Printer configuration: Resolution above 600dpi Others: Bluetooth 4.0	USB2.0
Network Requirement	Ethernet (not required) Ethernet is necessary to realize FTP / DICOM / WORKLIST / HL7 function	Ethernet (not required) Ethernet is necessary to realize FTP / DICOM / WORKLIST / HL7 function

### Instructions for use

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments and so on.

**Warning:** Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

**Warning:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**Warning:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

**Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment (detail model name), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**If any:** a list of all cables and maximum lengths of cables (if applicable), transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES may be specified either generically (e.g. shielded cable, load impedance) or specifically (e.g. by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).

**If any:** the performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term "ESSENTIAL PERFORMANCE" need not be used).

### Technical description

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
2. Guidance and manufacturer's declaration - electromagnetic emissions and Immunity

Table 1

<b>Guidance and manufacturer's declaration - electromagnetic emissions</b>	
<b>Emissions test</b>	<b>Compliance</b>
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class A
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Applied

**EN** Table 2

<b>Guidance and manufacturer's declaration - electromagnetic Immunity</b>		
<b>Immunity Test</b>	<b>IEC 60601-1-2 Test level</b>	<b>Compliance level</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV Contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV Contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Power supply lines: ±2 kV
Surge IEC 61000-4-5	Line(s) to line(s) : ±1 kV. Line(s) to earth : ±2 kV. 100 kHz repetition frequency	Line(s) to line(s) : ±1 kV. 100 kHz repetition frequency
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0,5 cycle At 0°, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° et 315 ° 0% 1 cycle and 70% 25/30 cycles Single phase: at 0 0% 300 cycle	0% 0,5 cycle At 0°, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° et 315 ° 0% 1 cycle and 70% 25/30 cycles Single phase: at 0 0% 300 cycle
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150 KHz to 80 MHz : 3Vrms 6Vrms (in ISM bands) 80% Am à 1kHz
Radiated RF IEC61000-4-3	10 V/m 80MHz- 2,7GHz 80% AM to 1kHz	10 V/m 80MHz- 2,7GHz 80% AM to 1kHz
<b>NOTE:</b> UT is the a.c. mains voltage prior to application of the test level.		

Table 3

<b>Guidance and manufacturer's declaration - electromagnetic Immunity</b>						
Radiated RF IEC61000-4-3						
(Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)						
<b>Test Frequency (MHz)</b>	<b>Band (MHz)</b>	<b>Service</b>	<b>Modulation</b>	<b>Modulation (W)</b>	<b>Distance (m)</b>	<b>IMMUNITY TEST LEVEL (V/m)</b>
385	380 -390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	380 -390	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	Bande LTE 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, Bande LTE 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						

Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
1720	1 700-1 990	GSM 1800 ; CDMA 1900 ; GSM 1900 ; DECT ; Bande LTE 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth WLAN, 802,11 b/g/n, RFID 2450, Bande LTE 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802, 11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5240						
5785						

### Appendix G: Manufacturer Information

**Manufacturer:** Shenzhen ECGMAC Medical Electronics Co., Ltd

**Factory Address:**

2nd Floor of Block 2, Haoye Industrial Park, Tiegang Road, Xixiang Street, Baoan District, 518102 Shenzhen, China

**Tel:** +86 755-27697821/ 27697823/ 27948579

**Fax:** +86 755-27697823-616

**Website:** www.ecgmac.com

**Email:** info@ecgmac.com

### Appendix H: EU Representative Information

**Company Name :** Wellkang Ltd


**Company Address :**

Enterprise Hub, NW Business Complex, 1 Beraghmore Rd., Derry, BT48 8SE, Northern Ireland, UK

**Tél. :** +44(20)3287 6300, +44(33)3303 1126

**Fax :** +44(20)76811874

**Web :** www.wellkang.ltd.uk ; www.CE-marking.eu

# Spengler

Ref.:NU\_CARDIOMATE\_PC\_VA\_140921