
## Röhrenhalter S U185001

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03/08 ALF

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
### Instruction sheet

03/08 ALF

- 1 Base plate
- 2 Slot for accommodating the Helmholtz coils
- 3 Ramp for accommodating the auxiliary coil
- 4 LED
- 5 Terminal panel (backside)
- 6 Neck brace
- 7 Tail stock with plug connectors

### 1. Safety instructions

When the tube is in operation, the terminals of the tube may be at high voltages with which it is dangerous to come into contact.

- Only use safety experiment leads to connect up the circuitry.
- Only change circuits with power supply equipment switched off.
- Only exchange tubes with power supply equipment switched off.

### 2. Description

The tube holder accommodates the entire range of tubes U18550-U18571 designed for practical experiments, along with their accessories such as Helmholtz coils U185051 and an auxiliary coil U19106, and can also serve as the base for the discharge tube U18580.

The tube holder consists of a base plate with a holder for the tube set at an angle with terminal panel attached. The tail-stock for the neck of the tube contains an octal plug connector and a central bore. Five sockets are electrically connected to the terminal panel with its 4-mm safety sockets. A filament protection circuit is also integrated into

the neck brace to prevent excess voltage being applied to the cathode heater circuit and thus protecting the sensitive heater filament from damage. The relay cuts out at a filament voltage of about 10.5 V DC or about 8.5 V AC. It switches back on again when the filament voltage is once again below the limit. Excess voltage is indicated by a red LED. The base also contains a slot with a metric scale for accommodating Helmholtz coils either transversely to the tube and up to 150 mm away from it or in standard Helmholtz configuration (marked H). The ramp at the front of the tube holder accommodates the Helmholtz coils in an axial configuration and also provides a base for the auxiliary coil. The tube holder is mounted on three rubber feet and is thus unable to slide.

Terminal	Tube pin	Function	Comment
A1	1	Anode	Connected directly
F3	3	Cathode heater	Filament protection circuit
F4	4	Cathode heater	Filament protection circuit connected to C5 via a 33 kΩ resistor
C5	5	Depends on the tube e.g. cathode	Connected directly and to F4 via a 33 kΩ resistor
G7	7	Depends on the tube e.g. grid	Connected directly

### 3. Technical data

Dimensions: 130 x 190 x 250 mm<sup>3</sup> approx.  
 Weight: 0.570 kg approx.

## 4. Operation

### 4.1 Setting up and removing tubes

- Press tube gently into the stock and push until the pins are fully inserted. Take note of the unique positions of the guide pin (refer to Fig. 1).
- To remove the tube, apply pressure with the index finger of the right hand on the guide pin until the pins loosen, then pull out the tube.

### 4.2 Setting up the Helmholtz coils for a transverse field

- Insert the coils into the middle of the coil slot and push them out as far as they go. Make sure the connectors point outwards (refer to Fig. 2).
- Insert the hot cathode tube into the holder as described above.
- For Helmholtz configuration  $d = r$  the outer edges of the coil base should be flush with the dashed lines.

### 4.3 Setting up the Helmholtz coils for an axial field

- Insert the hot cathode tube into the holder as described above.
- Slot the base of the coil into the groove from the front making sure that the connectors point forwards.
- When using both coils, attach the base of the second coil to the plug of the first coil (refer to Fig. 3a).

### 4.4 Setting up the auxiliary coil

- Connect up the coil using experiment leads (refer to Fig. 3b).
- Place the coil onto the ramp of the tube holder so that the plugs slide into the corresponding slot.
- The leads may be trailed through the mouth of the ramp.
- Insert the hot cathode tube into the holder as described above.


Fig. 1 Setting up a tube


Fig. 2 Setting up the coils for a transverse field


Fig. 3a Setting up a coil for an axial field


Fig. 3b Setting up the auxiliary coil

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