

SHK 18



SHK 18 Coussin Chauffant

Arrêt automatique au bout d'env. 90 minutes
Régulation électronique précise de la température
3 réglages de température
Ecran éclairé
Protection contre la surchauffe
Cordon détachable
Housse en coton amovible et lavable à la machine
Öko-Tex 100
Les textiles utilisés pour ce produit répondent aux exigences écologiques et rigoureuses de la certification Oeko-Tex Standard 100, et vérifié par l'Institut de recherches Hohenstein.
100 watts
Dimensions: 40 x 30 cm
Garantie 2 ans
UE: 5 / Colisage: -
EAN-Nr.: 4211125 24501 6
Art.-Nr.: 245.01

SHK 18 Heating pad

Automatic switch off after approx. 90 min.
Electronic temperature control
3 temperature settings
Illuminated display
Overheat protection
Detachable cord
Removeable cotton cover - washable
Oeko-Tex 100
The textiles used for this device meet the stringent human ecological requirements of Oeko-Tex Standard 100, as verified by Hohenstein Research Institute.
100 watt
Product measurements: 40 x 30 cm
2 year guarantee
Sales unit: 5 / Shipping carton: -
EAN no.: 4211125 24501 6
Item no.: 245.01