

# TERi™ Androgynous Geriatric Trainer

Product Name	<b>TERi™ Androgynous Geriatric Trainer</b>
Addresses Training Needs/Curricula	General capability carefully selected to support practice on: <ul style="list-style-type: none"> <li>- Essentials of Patient Care</li> <li>- Movement Assistance and Fall Prevention</li> <li>- Grooming and Daily Living Assistance</li> <li>- Basic and Advanced Nursing Skills</li> <li>- Cardiovascular and Respiratory Emergencies</li> </ul> Suitable for long-term facilities and emergency room training.
Configuration	<ul style="list-style-type: none"> <li>• Light or medium skin tones.</li> <li>• Configurable to fit unique training needs and budgets.</li> <li>• A range of accessories available that enhances the baseline version maximizing versatility, training benefit and return on investment.</li> </ul>
Accessories	IV Arm - LF04303 BP Arm - LF04304
Consumables	<ul style="list-style-type: none"> <li>• IV Arm <ul style="list-style-type: none"> <li>- Skins</li> <li>- Veins Kit</li> <li>- IM Pads</li> <li>- Fluid Supply Bag</li> </ul> </li> <li>• IM Pads for Thighs and Buttocks</li> <li>• Urinary Catheterization: <ul style="list-style-type: none"> <li>- Male Genitalia Foreskins</li> <li>- Foley Catheter</li> <li>- Urinary Catheterization Bladders</li> </ul> </li> <li>• Lung Bags</li> <li>• Stomach Bags</li> <li>• Stoma Bags</li> <li>• Gastrostomy Bags</li> <li>• Simulated Blood</li> <li>• Lubricant</li> </ul>
Product Code/SKU	Geriatric Patient Care Trainer (Light) - LF04300 Geriatric Patient Skills Trainer (Light) - LF04301 Geriatric Patient Simulator (Light) - LF04302
Dimensions/Weight	Dimensions: 5" feet Ship Weight: 65lbs



Be  
**READY**



**Nasco Healthcare Inc.**  
 16 Simulaids Drive, Saugerties, NY, 12477 USA  
 Phone: 1-833-NASCOHC (627-2642) Email: info@nascohealthcare.com  
 www.nascohealthcare.com



# Be Ready when called into action.

## TERi™ Androgynous Geriatric Trainer

A comprehensive male/female elderly patient care trainer for medical simulation. Evolving from our GERi™/KERi™ full-body manikins, it has a more realistic look and feel accurately representing the human anatomy. Weight is distributed to represent a real patient for lifting and carrying.

Each version is affordable and durable with upgradable items.

- **Geriatric Patient Care Trainer<sup>1</sup>** – General patient care, daily living assistance simulation platform.
- **Geriatric Patient Skills Trainer<sup>2</sup>** – Physical skills practice simulation platform with injection/IV, catheterization, cardiovascular and respiratory support.
- **Geriatric Patient Simulator<sup>3</sup>** – Comprehensive and complete simulation platform including patient monitoring capability.



## TERi™ Androgynous Geriatric Trainer



**PROM: Passive Range of Motion<sup>1, 2, 3</sup>**  
Seating, Laying, Finger/Toe Manipulation with No Effort from Patient

**Cardiovascular Support<sup>2, 3</sup>**  
Carotid Pulses, CPR with Quality Performance Metrics & Auscultation



**Intravenous, Injection, Subcutaneous Injection Sites<sup>2, 3, optional on 1</sup>**  
Multiple Locations Equipped with Standard Arm IM Injection Pads




**Nasogastric Tube Procedures<sup>2, 3</sup>**  
Placement, Lavage, Gavage, Feeding and Suctioning



**Enhanced Aesthetics<sup>1, 2, 3</sup>**  
Soft Silicone Skin with Noticeable Normal and Cancerous Skin Moles, Stage 1 Sacral Ulcer and Reddened Skin Folds

**Respiratory Support<sup>2, 3</sup>**  
Manual Ventilation with Chest Rise, Oxygen Administration, Tracheostomy Care, Airway/Trach Suctioning, Airway Intubation, Ventilation and Auscultation



**Grooming & Daily Living Assistance<sup>1, 2, 3</sup>**  
Bathing, Clothing Changes, Denture Placement/Removal, Hearing Aid Placement/Removal, Eye Irrigation & Optic Drops, Ostomy & Urinary Catheter Care, Hair Care



**Indwelling Male & Female Urinary Catheterization<sup>2, 3</sup>**  
Includes Support for Urethral and Suprapubic Catheters

**Advanced Skills<sup>2, 3</sup>**  
Urinary Catheterization, Enema Administration, Pelvic and Prostate Examination

This manikin is intended to be used as an instructional aide for teaching purposes. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies, or function as needed.